



Jul 12, 2019

## SPARC Serotonin 2A Receptor (5-HT2AR) Immunohistochemistry Protocol in Rat Tissues Labeled with Cholera Toxin B-fragment

Elisa Gonzalez-Rothi<sup>1</sup>, Yasin Seven<sup>1</sup>, Latoya Allen<sup>1</sup>, Marissa Ciesla<sup>1</sup>, Gordon Mitchell<sup>1</sup>

<sup>1</sup>University of Florida

1 Works for me dx.doi.org/10.17504/protocols.io.2khgct6



Elisa Gonzalez-Rothi

### ABSTRACT

This protocol describes the immunofluorescent labeling technique used to identify serotonin 2A receptor expression in CtB-labelled phrenic motor neurons and within a defined region of interest surrounding phrenic motor neurons.

- 1 Day 1: primary antibodies required:  
5-HT2AR: Rabbit anti-5-HT2AR (Immunostar #24288)  
Cholera toxin B-fragment: Goat anti-CT-B (Millipore #227040)
- 2 Place 40um transverse spinal cord sections into 1xPBS in 12 well cell culture plates
- 3 5x washes in 1xPBS for 5 minutes each at room temperature
- 4 Blocking: place tissues in 5% Normal Horse Serum (NHS) in 1xPBS-Triton (0.1%) +0.5% bovine serum albumin (BSA) for 60 minutes at room temperature
- 5 Primary Antibody Incubation: Incubate tissues in: 5%NHS in 1xPBS-Triton (0.1%), Rabbit anti-5-HT2AR (1:1000), and Goat anti-Ct-B (1:2500) for 48 hours at 4 degrees C
- 6 Day 2: secondary antibodies required:  
AlexaFluor 488: donkey anti-goat (Invitrogen, Ref#A11055)  
AlexaFluor 594: donkey anti-rabbit (Invitrogen, Ref#A21207)
- 7 5x washes in 1xPBS for 5 minutes each at room temperature
- 8 Secondary Antibody Incubation: Incubate tissues in: 5%NHS in 1xPBS-Triton (0.1%), donkey anti-goat (1:1000), and donkey anti-rabbit (1:1000) for 2 hours at room temperature
- 9 5x washes in 1xPBS for 5 minutes each at room temperature
- 10 Mount tissues on Superfrost Plus microscope slides
- 11 Allow slides to dry overnight

## 12 Coverslip with VectaShield Antifade Hard Set Mounting Medium (Cat#:H-1400)



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited