

PBS Buffer Solution

Dr. Steven Wilhelm

Abstract

Please contact Dr. Steven Wilhelm (wilhelm@utk.edu) for additional information regarding this protocol.

Citation: Dr. Steven Wilhelm PBS Buffer Solution. **protocols.io**

dx.doi.org/10.17504/protocols.io.in9cdh6

Published: 14 Jul 2017

Protocol

Step 1.

Add dH₂O to a clean, acid washed bottle

 **AMOUNT**

80 ml Additional info:

Step 2.

Add NaCl

 **AMOUNT**

8 g Additional info:

 **REAGENTS**

 Sodium chloride [View](#) by [P212121](#)

Step 3.

Add KCl

 **AMOUNT**

0.2 g Additional info:

 **REAGENTS**

 Potassium chloride [View](#) by [P212121](#)

Step 4.

Add Na₂HPO₄

 **AMOUNT**

1.44 g Additional info:

 **REAGENTS**

Sodium Phosphate dibasic S373-500 by [Fisher Scientific](#)

Step 5.

Add KH_2PO_4

AMOUNT

0.24 g Additional info:

REAGENTS

✓ Monopotassium phosphate by Contributed by users

Step 6.

Adjust the pH to 7.4 using HCl

Step 7.

Add dH_2O to a total volume of 100 mL