

HiFi DNA Assembly (NEB)

Josh Timmons

Abstract

This is the protocol for DNA Assembly using the NEBuilder® HiFi DNA Assembly Master Mix (E2621).

Citation: Josh Timmons HiFi DNA Assembly (NEB). [protocols.io](https://doi.org/10.17504/protocols.io.dn85hv)

[dx.doi.org/10.17504/protocols.io.dn85hv](https://doi.org/10.17504/protocols.io.dn85hv)

Published: 22 Aug 2015

Guidelines

Optimal Quantities

NEB recommends a total of 0.03–0.2 pmols of DNA fragments when 1 or 2 fragments are being assembled into a vector, and 0.2–0.5 pmols of DNA fragments when 4–6 fragments are being assembled. Efficiency of assembly decreases as the number or length of fragments increases. To calculate the number of pmols of each fragment for optimal assembly, based on fragment length and weight, we recommend the following formula, or using the tool, [NEBcalculator](#).

$\text{pmols} = (\text{weight in ng}) \times 1,000 / (\text{base pairs} \times 650 \text{ daltons})$

50 ng of 5000 bp dsDNA is about 0.015 pmols

50 ng of 500 bp dsDNA is about 0.15 pmols

The mass of each fragment can be measured using the NanoDrop instrument, absorbance at 260 nm or estimated from agarose gel electrophoresis followed by ethidium bromide staining.

Materials

NEBuilder HiFi DNA Assembly Master Mix - 10 rxns [E2621S](#) by [New England Biolabs](#)

Protocol

Step 1.

Set up the following reaction on ice (to 20µl total volume):

Recommended DNA Ratio	Recommended Amount of Fragments Used for Assembly		
	2–3 Fragment Assembly*	4–6 Fragment Assembly**	Positive Control†
	vector:insert = 1:2	vector:insert = 1:1	
Total Amount of Fragments	0.03–0.2 pmols*	0.2–0.5 pmols**	10 µl
	X µl	X µl	

NEBuilder HiFi DNA Assembly Master Mix	10 μ l	10 μ l	10 μ l
Deionized H2O	10-X μ l	10-X μ l	0
Total Volume	20 μ l † †	20 μ l † †	20 μ l

PROTOCOL

E2621 DNA Assembly Reaction

CONTACT: [New England Biolabs](#)

NOTES

Josh Timmons 22 Aug 2015

NEU iGEM found that a HiFi DNA Assembly with 7 pieces, each at 0.05pmol, was successful. Plasmid fragments, post-column purification, were diluted to 0.1pmol/uL in water. 0.5uL of each fragment was then pipette into the reaction tube on ice.

Step 1.1.

Vector DNA

ANNOTATIONS

Shine Sun 19 Apr 2016

pLenti Guide puro

Weijun Liu 15 Aug 2016

pLentiCRISPRV2

Step 1.2.

Insert fragments DNA

Step 1.3.

NEBuilder HiFi DNA Assembly Master Mix

REAGENTS

 NEBuilder HiFi DNA Assembly Master Mix - 10 rxns [E2621S](#) by [New England Biolabs](#)

Step 1.4.

Deionized H2O

Step 2.

Incubate samples in a thermocycler at 50°C for 15 minutes (when 2 or 3 fragments are being assembled) or 60 minutes (when 4-6 fragments are being assembled).

DURATION

01:00:00

Step 3.

Following incubation, store samples on ice or at -20°C for subsequent transformation.

Step 4.

Transform NEB 5-alpha Competent E. coli cells (provided in the cloning kit or purchased separately from NEB) with 2 μ l of the assembled product, following the [transformation protocol](#).