



Oct 09, 2018

Working

# Dissociation of Nuclear Suspensions from Human Breast Tissues

Emi Sei<sup>1</sup>, Shanshan Bai<sup>1</sup>, Nicholas Navin<sup>1</sup><sup>1</sup>University of Texas MD Anderson Cancer Center

dx.doi.org/10.17504/protocols.io.t3aeqie

Nicholas Navin

## ABSTRACT

This protocol can be used to prepare nuclear suspensions from fresh or frozen tissue for single cell DNA or RNA sequencing experiments. For RNA analysis it is recommended to use fresh tissue samples, due to RNA degradation that occurs during freeze-thaw cycles, however for DNA analysis either frozen or fresh tissue can be used. The protocol has been tested on breast tissues, and many frozen tumor tissues including breast, colon prostate and pancreatic samples.

## TAGS

breast cancer



NL\_Breasttissue nuclei  
dissociation\_protocol\_  
NN.docx

## PROTOCOL STATUS

### Working

We use this protocol in our group and it is working

## GUIDELINES

This protocol is for fresh or frozen tissue.

## MATERIALS

NAME	CATALOG #	VENDOR
Flowmi™ Cell Strainer 40 µm	H13680-0040	Bel-Art
DNA LoBind Tube 1.5ml	022431021	Eppendorf
Ultrapure BSA	AM2616	Ambion
RNase Inhibitor (40U/µl)	N2615	Promega

## MATERIALS TEXT

### Consumables:

- Scalpels, 1.5ml tubes, 50ml conical, syringe, PBS.

## SAFETY WARNINGS

Please refer to the SDS (Safety Data Sheet) for safety warnings and hazard information.

## BEFORE STARTING

- NST-DAPI: Mix 800 ml of NST solution (146 mM NaCl, 10 mM Tris base (pH 7.8), 1 mM CaCl<sub>2</sub>, 21 mM MgCl<sub>2</sub>, 0.05% (wt/vol) BSA and 0.2% (vol/vol) Nonidet P-40) with 200 ml of DAPI solution (106 mM MgCl<sub>2</sub> and 10 mg of DAPI). Filter-sterilize the

solution and store it at 4°C in the dark for up to 1 year.

- Prepare chilled (4°C) Nuclei Wash and Resuspension Buffer: 1X PBS with 1.0% BSA and 0.2U/μl RNase Inhibitor.

### Tissue preparation

- 1 Cut fresh or frozen tissue (about 1x1x1cm) and mince it with a scalpel in a 10-cm Petri dish with **2 ml NST-DAPI buffer** for **00:10:00** (10-15 min) until tissue chunks are no longer visible.
- 2 Filter the mixture through a 40 μl flowmi into 1.5 ml DNA LoBind tubes. Place the tube on ice.

### Centrifugation

- 3 Centrifuge nuclei at **4 °C** 300-500 xg for **00:05:00**.

### Wash

- 4 Remove the supernatant and resuspend with **1 ml cold Nuclei Wash and Resuspension Buffer**. (1/3)
- 5 Wash again with **1 ml cold Nuclei Wash and Resuspension Buffer**. (2/3)
- 6 Wash again with **1 ml cold Nuclei Wash and Resuspension Buffer**. (3/3)

### Counting nuclei

- 7 Resuspend pellet in cold Nuclei Wash and Resuspension Buffer. Count nuclei (Make a 1:1 solution of Hoechst and PI, pipette **2 μl** of this mix to **10 μl filtered nuclei** and count in the Countess II FL).
- 8 Adjust concentration to 700-1200 nuclei/μl.
- 9 Proceed immediately with 10X genomics RNA experiments or single cell DNA sequencing experiments.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited