

Bacteriophage Cryopreservation

Dr. Steven Wilhelm

Abstract

Please contact Dr. Steven Wilhelm (wilhelm@utk.edu) for additional information regarding this protocol.

Citation: Dr. Steven Wilhelm Bacteriophage Cryopreservation. **protocols.io**

dx.doi.org/10.17504/protocols.io.iyyccfw

Published: 24 Jun 2017

Protocol

Preparation

Step 1.

Vortex a freshly lysed culture for 10 sec

Preparation

Step 2.

Add 1 mL of phage into each cryogenic vial

 **AMOUNT**

1 ml Additional info:

Preparation

Step 3.

Add 1 mL of 20% glycerol to each vial with phage

 **AMOUNT**

1 ml Additional info:

 **REAGENTS**

Glycerol MRGE-4002 by growcells.com

Preparation

Step 4.

Invert tubes several times

Preparation

Step 5.

Place immediately in a -80°C freezer

Transferring frozen virus particles

Step 6.

Remove cryogenic vials one at a time from the freezer

Transferring frozen virus particles

Step 7.

Using a sterile pipette tip, scrape a small amount of frozen virus particles into a 20 mL mid-log culture

Transferring frozen virus particles

Step 8.

Place cryogenic vials immediately back into the freezer

Transferring frozen virus particles

Step 9.

Allow to incubate and lyse culture under constant light conditions appropriate for your specific culture