



Sep 22, 2019

## Plaque assay

Marijn Ceelen<sup>1</sup>

<sup>1</sup>iGEM Wageningen

1 Works for me dx.doi.org/10.17504/protocols.io.7i8hkhw

iGEM Wageningen 2019

Marijn Ceelen  
Wageningen University

### ABSTRACT

This is a protocol for the quantification of phage lambda and T7 titers by counting plaques.

### MATERIALS

NAME	CATALOG #	VENDOR
Liquid LB medium		
LB agar	<a href="#">View</a>	

- 1 Streak LB plate with E. coli (strain should be susceptible to phage lambda. Strains that include a prophage lambda are likely to be resistant. Strains LE392 and DH10B have been used with this protocol) and incubate overnight at  $37^{\circ}\text{C}$ .
- 2 Pick a colony from this plate and use it to inoculate  $10\text{ ml}$  of LB. Incubate this culture at  $37^{\circ}\text{C}$  until OD reaches 2-3. Overnight culture is recommended.
- 3 Make a range of dilutions. Mix  $100\text{ }\mu\text{l}$  of phage dilution with  $200\text{ }\mu\text{l}$  of E. coli culture. Then incubate for 10 minutes at room temperature. After  $00:10:00$  add  $3\text{ ml}$  of liquid soft LB agar (LB agar with 0.7% agar) ( $50^{\circ}\text{C}$ ).
- 4 Pour the resulting mixture on LB plates (preheat plates at  $37^{\circ}\text{C}$ ). Spread the mixture on the plate by moving it. Make sure to work quickly to avoid clumps of solidified agar.
- 5 Incubate for  $00:15:00$  at room temperature. Afterwards turn the plates over and incubate overnight at  $37^{\circ}\text{C}$ .
- 6 Use plates with 30-300 plaques to determine phage concentration. Calculate the Plaque Forming Units (PFU)/mL by the following formula:  $\text{PFU/mL} = N \times 1/DF \times 1/V$ .  
N is the number of plaques of lysis counted on the plate (expressed as PFU); DF is the dilution factor and V is the volume of phage dilution poured on the plate.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited