



## Platereader workflow with *V. natriegens*

Daniel Stukenberg<sup>1</sup>

<sup>1</sup>Philipps- university Marburg

[dx.doi.org/10.17504/protocols.io.uy9exz6](https://dx.doi.org/10.17504/protocols.io.uy9exz6)

iGEM Marburg 2018

 Daniel Stukenberg 

### ABSTRACT

This protocol provides a workflow for platereader measurements with *V. natriegens*.

### TAGS

plate reader

GFP

Show tags

### PROTOCOL STATUS

#### Working

We use this protocol in our group and it is working

### MATERIALS TEXT

transparent and black 96-well plate

#### Sample preparation

- 1 Aliquot 50  $\mu$ L LBv2 in 1.5 mL reaction tubes
- 2 Transfer material from glycerol stock into these reaction tubes
- 3 Prepare transparent 96 well plat with 190  $\mu$ L LBv2
- 4 Use 10  $\mu$ L of LBv2 with cells from glycerol stock to inoculate the 96 well plate. (Carry out experiment with four technical replicates)

#### Preculture

- 5 Incubate 96 well plate in a platereader (Protocol: Preculture)

Preculture:

37 cycles

600 s cycle time

Shaking: Double orbital, 500 rpm

Protocols: OD + optional protocols

OD600 (Absorbance):  
Wavelength: 600 nm  
Settling Time: 0.5 s  
No. of flashes: 30

6 Dilute grown cultures 1:40 (195  $\mu$ L LBv2 + 5  $\mu$ L culture) in a sterile 96 well plate

7 Apply Shaking protocol

Shaking (Absorbance):  
Shaking: Double Orbital, 300s, 700 rpm  
Measure OD600

8 Dilute grown cultures 1:50 (196  $\mu$ L LBv2 + 4  $\mu$ L culture) in a black sterile 96 well plate

#### Measurement

9 Start protocol for measurement in the plate reader (Protocol: OD+Lux, OD+GFP, OD+RFP)

OD + Lux  
121 cycles  
300 s cycle time  
Shaking: Double orbital, 500 rpm  
Protocols: OD + luminescence

OD600 (Absorbance):  
Wavelength: 600 nm  
Settling Time: 0.5 s  
No. of flashes: 30

Luminescence:  
Gain: 4000  
Measurement interval time: 1 s



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited