



Jul 24, 2019

Prepare NGM no peptone plates V.2

Cancer Research UK / Wellcome Gurdon Institute media kitchen¹¹Cancer Research UK / Wellcome Trust Gurdon Institute

1

Works for me

dx.doi.org/10.17504/protocols.io.5sng6de

Cristian Riccio



ABSTRACT

Prepare NGM plates without peptone.

Peptone is necessary for bacterial growth and thus by not adding peptone, one can control for the number of bacteria present on an NGM plate by seeding a known number of bacteria.

MATERIALS

| NAME ▾ | CATALOG # ▾ | VENDOR ▾ |
|--|---------------------------|-----------------------------------|
| double distilled water (ddH2O) | | |
| Sodium chloride meets analytical specification of Ph.Eur Fisher Chemical | S/3160/65 | Fisher Scientific |
| SYCHEM autoclave | View | Syschem |
| Agar | AGA02 | Formedium |

SAFETY WARNINGS

Make sure you know how to use the autoclave before starting this protocol.

1 [NGM No Peptone agar_5L.xls](#)

2

| Ingredients | Quantity | |
|----------------------|------------------|------|
| NaCl | | 15g |
| Double distilled H2O | 4.8L | |
| Agar | per 1L bottle | 17g |
| Agar | per 500ml bottle | 8.5g |

3

| |
|--|
| Measure approx 4.8L double distilled H2O in 5L bell jar with a magnetic flea. |
| Add NaCl |
| Stir until all solutes are dissolved. |
| Dispense approx 972ml NGM media per 1L bottle OR 486ml NGM into a 500ml bottle. |
| Dispense 17g of agar and a magnetic flea to each 1L bottle OR 8.5g of agar and a magnetic flea to each 500ml bottle. |
| Label, date and autoclave. |



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited