

# Euplotes focardii RNA extraction

Angela Piersanti

## Abstract

**Citation:** Angela Piersanti Euplotes focardii RNA extraction. **protocols.io**

dx.doi.org/10.17504/protocols.io.hhwb37e

**Published:** 01 Apr 2017

## Protocol

### Step 1.

Filter and collect Euplotes focardii cells (ca.40000 cells) by centrifuging at 4°C for 5 min at 3000 rpm.

### Step 2.

Dissolve the cell pellet in 750 µl Trizol reagent solution (3 volumes of Trizol to 1 volume of cells), vortex for 2 min.

### Step 3.

Incubate 5 min at room temperature, add 200 µl of chloroform to the solution and vortex for 15 sec.

### Step 4.

Incubate 15 min at room temperature then centrifuge the sample at 4°C for 20 min at 12000 rpm.

### Step 5.

Move the supernatant containing RNA in a new tube and precipitate the RNA by adding isopropanol vol1:1.

### Step 6.

Incubate 10 min at room temperature, collect the RNA by centrifugation at 4°C for 30 min at 13000 rpm and discard the supernatant.

### Step 7.

Wash twice the RNA pellet with 1 ml of ethanol 75% (in DPEC water) by centrifuging 5 min at 8000 rpm and then let the pellet dry.

### Step 8.

Resuspend the RNA pellet in 20-30 µl of DPEC water and incubate 10 min at 55-60°C.

**Step 9.**

Treat the RNA sample with DNase I to remove genomic DNA.

**Step 10.**

Treat the RNA sample with DNase I to remove genomic DNA.