

Colony PCR on *V. natriegens* Version 2

Daniel Bauersachs

Abstract

Citation: Daniel Bauersachs Colony PCR on *V. natriegens*. **protocols.io**

dx.doi.org/10.17504/protocols.io.pundnve

Published: 30 Apr 2018

Protocol

Step 1.

5 µl	reaction buffer
0,5 µl	dNTPs
0,5 µl	primer forward
0,5 µl	primer reverse
0,1 µl	one-taq polymerase
18,4 µl	A.dest.
= 25 µl	

Now you have to pick one single colony with a pipette tip and put it into your reaction mix. After 15 seconds you can take the tip out and start the PCR.