



May 17,  
2019

Working

## Yale - Blood Urea Nitrogen [↗](#)

John Stack<sup>1</sup>, Gary Cline<sup>1</sup>

<sup>1</sup>Yale University

[dx.doi.org/10.17504/protocols.io.y3qfymw](https://doi.org/10.17504/protocols.io.y3qfymw)

Mouse Metabolic Phenotyping Centers  
Tech. support email: [info@mmpc.org](mailto:info@mmpc.org)

Lili Liang

### ABSTRACT

#### Summary:

Procedure used to measure the concentration of Blood Urea Nitrogen(BUN) in blood, plasma, and serum. Urea is determined by the enzymatically coupled reactions of urease (to form ammonia) and glutamate dehydrogenase (conversion of ammonia and glutamate to glutamine with oxidation of NADH to NAD). The rate of NAD formation is monitored by the change in absorbance at 340 nm.

### EXTERNAL LINK

<https://mmpc.org/shared/document.aspx?id=203&docType=Protocol>

### MATERIALS

NAME	CATALOG #	VENDOR
BUN liquid Reagent	R84533	Prolabs(cliniqa)
Multi Analyte Calibrator	R60010	Prolabs(cliniqa)
Assayed Control Serum 1	R83082	Prolabs(cliniqa)
Assayed Control Serum 2	R83083	Prolabs(cliniqa)

### MATERIALS TEXT

#### Reagent Preparation:

**BUN liquid Reagent:** As supplied by vendor

**Multi Analyte Calibrator:** Add the appropriate amount of water (6.5mL) to the chemical control bottle. Invert to mix, allowing 15 minutes for the reagent to settle.

**Assayed Control Serum 1:** Add the appropriate amount of water (6.5mL) to the chemical control bottle. Invert to mix, allowing 15 minutes for the reagent to settle.

**Assayed Control Serum 2:** Add the appropriate amount of water (6.5mL) to the chemical control bottle. Invert to mix, allowing 15 minutes for the reagent to settle.

### BEFORE STARTING

*Analysis by automated system Cobas Mira Plus*

- 1 Calibrate Cobas for BUN analysis by running a multi analyte standard and two assayed control serums.

- 2 Sample handling as performed by the Cobas Mira Plus.
- a) Cobas pipettes 2  $\mu\text{L}$  of sample into a cuvette slot.
  - b) Absorbance is measured at 340 nm.
  - c) Add 200  $\mu\text{L}$  of BUN liquid reagent.
  - d) Mixture is incubated at 37°C for 10 minutes.
  - e) Absorbance is measured at 340 nm. Change in absorbance is calculated.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited