



Dec 03, 2019

Domestication of L0 parts for Loop type IIS (Bsal and SapI)
Forked from [Domestication of L0 parts for Loop type IIS \(Bsal and SapI\)](#)

[Eftychis Frangedakis](#)¹, [marta tomaselli](#)¹, [Marius Rebmann](#)¹, [Susana Sauret-Gueto](#)¹

¹Plant Sciences, University of Cambridge, OpenPlant

1 *Works for me* [dx.doi.org/10.17504/protocols.io.92jh8cn](https://doi.org/10.17504/protocols.io.92jh8cn)

OpenPlant Project

 **Susana Sauret-Gueto**
Plant Sciences, University of Cambridge, OpenPlant   


ABSTRACT

Domestication and cloning of L0 parts into pUAP4 for [Loop type IIS assembly](#)

MATERIALS

NAME	CATALOG #	VENDOR
Phusion high-fidelity PCR kit	F553S	Thermo Scientific
MinElute Gel Extraction Kit	28606	

1



Domestication of internal BsaI and SapI sites is necessary when a DNA sequence of interest contains a BsaI or SapI recognition site. For domestication, overlapping PCR can be used. Briefly, two PCR fragments are amplified upstream (fragment 1, with primers 1 and 2) and downstream (fragment 2, with primers 3 and 4) of the sequence to domesticate, which will be part of the reverse primer (primer 2) of fragment 1 and of the forward primer (primer 3) of fragment 2. Both primer 2 and 3 will be specially designed with a single nucleotide mismatch to alter the sequence to domesticate (taking care to not alter amino acid composition if the region is in the coding sequence). Primers will also contain SapI recognition sites to allow the two fragments to be ligated together into pUAP4 using SapI type IIS assembly. Primer 1 and primer 4 are designed according to the guidelines for [cloning of L0 parts into pUAP4](#). Blue arrows: SapI recognition site. Blue dashed lines: SapI cleavage site. Red arrows: BsaI recognition site. CmR: chloramphenicol bacterial resistance cassette. LacZ: lacZα cassette for blue-white screening of colonies.

PCR and amplified fragments purification

- 2 Amplify the parts using your preferred polymerase (e.g. Phusion, Thermo Scientific™) following the manufacturer instructions.
- 3 Run your PCR product on agarose gel.
- 4 Extract the amplified fragments using a PCR purification kit (e.g. Qiagen MinElute™).

- 5 Set up the assembly of the fragments as in [cloning of L0 parts into pUAP4](#)

[pUAP4]= plasmid length/200

[fragment]= fragment length/100

The formulas above will give you the required concentration in ng/μl.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited