



Jul 05, 2019

## Cryopreservation of Mammalian Cells (Suspension)

Andrew Crowley<sup>1</sup><sup>1</sup>Dartmouth College

Working

[dx.doi.org/10.17504/protocols.io.47pgzmn](https://doi.org/10.17504/protocols.io.47pgzmn)

Ackerman Lab

Andrew Crowley  
Dartmouth College, Ackerman Lab

### ABSTRACT

A protocol for the preservation of suspension-type cells by freezing.

The protocol has been successfully used on:

- THP-1 monocytes

### MATERIALS

NAME	CATALOG #	VENDOR
RPMI 1640 with L-glutamine	10040CV	Corning
MEM Nonessential Amino Acids	25025CI	Corning
Sodium Pyruvate	25000CI	Corning
Fetal Bovine Serum	S1620-50	BioWest

### BEFORE STARTING

Formulate and warm cryopreservation media consisting of:

- 50% fetal bovine serum
- 40% culture medium (may contribute additional FBS)
- 10% DMSO

- 1 Measure the density and viability of the cells

TC20 Automated Cell Counter  
automated cell counter

BioRad 1450102



- 2 Gently pellet the cells via centrifugation

250 x g 00:07:00

- 3 While centrifuging, use the **live** cell counts from **step 1** to calculate the total volume required to resuspend the cells @

[M]2500000 cells/mL


- 4 Gently aspirate the culture media from the pellet
- 5 Gently resuspend the cell pellet in cryopreservation media:
  - [M]50 Volume Percent fetal bovine serum
  - [M]40 Volume Percent cell culture media
  - [M]10 Volume Percent DMSO

Use the media volume calculated in **step 3**

- 6 Aliquot cells into labeled cryo vials

 1 ml

- 7 Place cryo vials into a **room temp** Mr. Frosty freezing container



Mr. Frosty Freezing Container  
freezing container

ThermoFisher Scientific 5100-0001 [↗](#)

- 8 Freeze the container overnight

 -80 °C  16:00:00 (ie. overnight)

- 9 Transfer the vials to liquid nitrogen for long term storage

- 10 Suggested cell culture media for **step 5**

---

step case

---

### THP-1s

RPMI-1640, supplemented with:

- sodium pyruvate (1 mM)
- non-essential amino acids (1x)
- fetal bovine serum (10%)

---

step case

---

### THP-1s

RPMI-1640, supplemented with:

- sodium pyruvate (1 mM)
- non-essential amino acids (1x)
- fetal bovine serum (10%)



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited