



Zoospore Isolation of *Aurantiochytrium limacinum*

Mariana Rius¹

¹State University of New York at Stony Brook

[dx.doi.org/10.17504/protocols.io.u73ezqn](https://doi.org/10.17504/protocols.io.u73ezqn)

Protist Research to Optimize Tools in Genetics (PROT-G)

Collier Lab

 Mariana Rius   

ABSTRACT

This protocol results in the isolation of single nucleated zoospores from a culture of *Aurantiochytrium limacinum*. The point is to recover single nucleus colonies rather than single cell colonies (which may contain varying nuclei - particularly following transformation).

TAGS

isolation

zoospores

PROTOCOL STATUS

Working

We use this protocol in our group and it is working

MATERIALS TEXT

Filtration Pump (or other way to pull a vacuum)

Filter holder (mftr#: Thermo Scientific Nalgene 300-4000)

5 um nylon filter (catalog#: R50SP04700)

15 ml culture tube

GPY media (+ zeocin 100 ug/ml)

- 0.5% Yeast Extract, 1% Peptone, 3% D+-Glucose, 1.8% instant ocean

GPY zeocin plate (+ zeocin 100 ug/ml + ampicilin 100 ug/ml)

Prepare Culture

- 1 Inoculate 500 ul GPY media (microcentrifuge tube) with toothpick scrape of a colony. If inoculating a zeocin resistant colony, add 5 ul zeocin (for 100 ug/ml). Vortex/ resuspend cells. Pipette 4.5 ml of GPY into a 15 ml glass culture tube. Add 500 ul inoculum. Incubate at 28 C, 170 rpm for 48 h, or until microscopy reveals the presence of zoospores.

Filter Culture

- 2 Filter culture using a 5 um nylon filter membrane, nalgene filter holder, and filtration pump. Transfer flowthrough to a microcentrifuge tube and centrifuge 5 min 11000 rpm.

Plate Zoospores

- 3 Resuspend cell pellet in 100 ul of ½ Artificial Seawater (ASW) or media. Make a 10X dilution by adding 10 ul of cells to 90 ul of ½ ASW or media. Spread onto 2% agar GPY (+ zeocin +amp) plate. Incubate at 28 C for 48 - 72 h until colonies appear.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited