

# Seamless editing of sgRNAs cloning mixture

Alexandre Paix

## Abstract

**Citation:** Alexandre Paix Seamless editing of sgRNAs cloning mixture. **protocols.io**

[dx.doi.org/10.17504/protocols.io.dxb7im](https://dx.doi.org/10.17504/protocols.io.dxb7im)

**Published:** 06 Oct 2015

## Protocol

### Step 1.

Q5 2X master mix (12.5ul)

### Step 2.

Forward primers 10uM (1.25ul)

### Step 3.

Reverse primers 10uM (1.25ul)

### Step 4.

pDD162 from a 1.5ml bacterial culture miniprep (0.5ul)

### Step 5.

H2O (9.5ul)