



Working

DNA extraction and Nested PCR [↗](#)

PLOS Neglected Tropical Diseases

Baiyan Gong¹ , Yaming Yang² , Xiaohua Liu¹, Jianping Cao^{3,4,5,6,7,8}, Meng Xu^{3,4,5,6,7,8}, Ning Xu^{3,4,5,6,7,8}, Fengkun Yang¹, Fangwei Wu², Benfu Li², Aiqin Liu^{*1}, Yujuan Shen^{*3,4,5,6,7,8}

¹Department of Parasitology, Harbin Medical University, Harbin, Heilongjiang, China, ²Yunnan Institute of Parasitic Diseases, Puer, Yunnan, China, ³National Institute of Parasitic Diseases, Chinese Center for Disease Control and Prevention, ⁴Chinese Center for Tropical Diseases Research, ⁵WHO Collaborating Centre for Tropical Diseases, ⁶National Center for International Research on Tropical Diseases, Ministry of Science and Technology, ⁷Key Laboratory of Parasite and Vector Biology, MOH, ⁸Shanghai, China

[dx.doi.org/10.17504/protocols.io.yvjfw4n](https://doi.org/10.17504/protocols.io.yvjfw4n)

 Baiyan Gong 

ABSTRACT

This is the protocol for the Nested PCR assay

EXTERNAL LINK

<https://doi.org/10.1371/journal.pntd.0007356>



  410bp.pdf

GUIDELINES

All primer sequences are listed below:

The first round of PCR was carried out using EbGeno-fe (5'-TTC AGA TGG TCA TAG GGA TG-3') as the forward primer and EbGeno-re (5'-ATT AGA GCA TTC CGT GAGG-3') as the reverse primer, which together amplify a 465-bp specific fragment.

A second round of PCR was performed using the forward primer EbGeno-fi (5'-TCG GCT CTG AAT ATC TAT GG-3') and the reverse primer EbGeno-ri (5'-ATT CTT TCG CGC TCG TC-3') in order to amplify a 410-bp internal fragment used for genotype specification.

MATERIALS TEXT

DNase/RNase-Free Deionized Water, primer, TaKaRa Taq DNA Polymerase (TaKaRa Bio Inc., Tokyo, Japan), 10×PCR Buffer, 10 mmol dNTP, 1.5 % agarose gel, GelStrain

- 1** DNA extraction
Genomic DNA of *E. bienersi* was extracted directly from 180 to 200 mg of each fecal specimen using the recommended procedures and the provided reagents by the manufacture of QIAamp DNA Stool Mini Kit (QIAGEN, Hilden, Germany). To obtain high yield of DNA, the lysis temperature was increased to 95 °C according to the manufacturer's suggestion. DNA was eluted in 200 µL of AE and stored at -20°C before it was used for PCR analysis.
- 2** PCR amplification
Amplification was accomplished in final volume of 25 µL, containing 17 µL of DNase/RNase-Free Deionized Water, 0.5 µL of each primer, 0.5 µL of TaKaRa Taq DNA Polymerase (TaKaRa Bio Inc., Tokyo, Japan), 2.5 µL of 10 × PCR Buffer, 2 µL of 10 mmol dNTP and 2 µL of DNA.
- 3** The target DNA undergoes the first run of PCR with the first set of primers.

 00:05:00 Initial denaturation

🔥 95 °C

4 ⌚ 00:00:40 Denaturation

🔥 94 °C

5 ⌚ 00:00:45 Annealing

🔥 53 °C

6 ⌚ 00:00:45 Extension/elongation

🔥 72 °C

7 [🔄 go to step #4 35 cycles](#)

8 ⌚ 00:04:00 Final extension

🔥 72 °C

9 The product from the first reaction undergoes a second run of PCR with the second set of primers.

⌚ 00:05:00 Initial denaturation

🔥 95 °C

10 ⌚ 00:00:35 Denaturation

🔥 94 °C

11 ⌚ 00:00:40 Annealing

🔥 55 °C

12 ⌚ 00:00:40 Extension/elongation

🔥 72 °C

13 [🔄 go to step #10 30 cycles](#)

14 ⌚ 00:05:00 Final extension

🔥 72 °C

15 Final hold

The final step cools the reaction chamber to 4°C for an indefinite time, and may be employed for short-term storage of the PCR products.

- 16 All secondary PCR products were subjected to electrophoresis in a 1.5 % agarose gel and visualized by staining the gel with GelStrain (TransGen Biotech., Beijing, China).



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited