

# CGAP MACS Live Dead Separation

Adam Hunter

## Abstract

**Citation:** Adam Hunter CGAP MACS Live Dead Separation. **protocols.io**

[dx.doi.org/10.17504/protocols.io.qz5dx86](https://dx.doi.org/10.17504/protocols.io.qz5dx86)

**Published:** 11 Jul 2018

## Protocol

### Step 1.

Material	Quantity	Supplier Info
15ml Falcon Tubes	3	Falcon (352097)
50ml Falcon Tubes	1	Falcon (352098)
MACS Dead Cell Removal Kit	1	Miltenyi Biotech (130-090-101)
Nuclease Free Water	19ml	Ambion (AM9939)
LS Columns	1	Miltenyi Biotech (130-042-401)
0.5ml DNA LoBind Eppendorf Tubes	1	Eppendorf (0030 108.035)
Trypan Blue	20ul	Fisher Scientific (11414815)
C-Chips	1	Cambridge Bioscience (DHC-N01-50)
PBS	10ml	GIBCO (14190-144)
Bovine Serum Albumin (BSA)	400ul	Sigma-Aldrich Co. Ltd (A7906-10G)

### Step 2.

A single-cell suspension should have been prepared previously and cells number and viability assessed using 1:1 trypan blue dilution.

- A viability percentage below 70-80% usually justifies using this Dead Cell Removal protocol.

### Step 3.

Remove required number of cells and place in a 15ml Falcon Tube.

- Required number of cells/total cells = volume required (ml).

### Step 4.

Prepare 20ml 1X Binding Buffer by adding 1ml 20X Binding Buffer Stock to 19ml Nuclease Free Water.

### Step 5.

Centrifuge cell suspension for 5min at 300g.

**Step 6.**

Remove supernatant.

---

**Step 7.**

Resuspend cell pellet in 100ul Dead Cell Removal MicroBeads per  $10^7$  cells.

---

**Step 8.**

Mix well and incubate for 15mins at room temperature.

---

**Step 9.**

When 5min of incubation remains, place MS column (if  $<2 \times 10^8$  cells) or an LS column (if  $<2 \times 10^9$  cells) on QuadroMACS Magnetic Cell Separator and run 500 $\mu$ l (MS column) or 3ml (LS column) 1X Binding Buffer through the LS column, using a waste 15ml Falcon Tube to catch the effluent.

---

**Step 10.**

When incubation is finished, add 1ml (MS column) or 3ml (LS column) 1X Binding Buffer to cells.

---

**Step 11.**

Run cell suspension through LS column on QuadroMACS Magnetic Cell Separator, using a 15ml Falcon Tube to catch effluent as the the live cell fraction.

---

**Step 12.**

When cells have passed through, run 4 x 500 $\mu$ l (MS column) or 4 x 3ml (LS column) 1X Binding Buffer through LS column on QuadroMACS Magnetic Cell Separator using the same falcon tube to catch effluent as the the live cell fraction.

---

**Step 13.**

Centrifuge cells at 500g for 5 min at 4°C. Resuspend in 0.5-1ml PBS + 0.04% BSA.

---

**Step 14.**

Count cells and viability using nucleocounter.

---

**Step 15.**

Resuspend in appropriate volume of 0.04% BSA in PBS to run in Chromium.

---