



Jun 11,
2019

IHC on fixed frozen cryo sections

Grace Burgin¹, Inbal Benhar¹

¹Broad Institute

Working

[dx.doi.org/10.17504/protocols.io.32egqbe](https://doi.org/10.17504/protocols.io.32egqbe)



Grace Burgin



GUIDELINES

All incubation steps in humidifier tray (slide box with PBS soaked paper towels)

All wash steps in coplin jars (no agitation)

- 1 30-60 min Protein Block
5% Normal Horse (or Donkey) Serum in PBS + 0.3% Triton-X
- 2 Overnight (~18-22hrs) Primaries at 4C in PBS + 0.3% Triton-X
- 3 Wash 2x5 min w/ PBS
- 4 2hr Secondaries at room temp in PBS + 0.3% Triton-X
- 5 Wash 2x5 min w/ PBS
- 6 Coverslip



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited