

# Isolation of leucocytes from human blood

Anna F. Zeitler, Luisa F. Jiménez-Soto

## Abstract

Protocol established by Anna F. Zeitler and published by Luisa F. Jiménez-Soto

This protocol is the final adaptation of protocols used in the laboratory of Prof. Rainer Haas. The following people contributed with their ideas and experience: Benjamin Busch, Bettina Vogl-Gebert and Luisa F. Jiménez-Soto

**Citation:** Anna F. Zeitler, Luisa F. Jiménez-Soto Isolation of leucocytes from human blood. **protocols.io**

dx.doi.org/10.17504/protocols.io.hjxb4pn

**Published:** 08 Jun 2017

## Protocol

### Extract human blood

#### Step 1.

Following ethical regulations and hygiene standards, whole blood was drawn from human volunteers into a collection tube containing at least 20u/ml of Heparin.

### Prepare blood for gradient centrifugation

#### Step 2.

- Dilute blood (12ml whole blood and 25 ml PBS/2mM EDTA, Mix by inversion).
- Add 35 ml of the diluted blood carefully onto 15 ml of Percoll™ Plus 55% isotonic (GE Healthcare Life Science). Avoid mixing cells and Percoll® Plus.



#### REAGENTS

Percoll™ Plus [17-5445-01](#) by [GE Healthcare](#)

### Centrifuge to separate the leucocytes from erythrocytes

#### Step 3.

Centrifuge the gradient in a swing rotor centrifuge at 400g for 20 minutes at 16°C without brake.

### Extract leucocytes from gradient

#### Step 4.

At the end of centrifugation, remove carefully the tube from the rotor. You should see the formation of three layers: i) Serum + Platelets; ii) Leucocytes, and iii) Erythrocytes.

Wash, collect and prepare leucocytes.

#### **Step 5.**

- Centrifuge the leucocyte solution in a swing-rotor centrifuge at 200 g for 5 min 16°C.
- Remove supernatant and repeat twice with PBS to remove any trace of Percoll Plus
- Resuspend cells in RPMI 1640 media complemented with 10% FBS.
- Count cells using Trypan Blue for viability count and seed on cell culture treated plates as needed ( $1 \times 10^6$  cells per well usually).

Incubate cells

#### **Step 6.**

Place cells at 37°C in a 5% CO<sub>2</sub> atmosphere. We used the cells in the next 24 hours to secure a close physiological behaviour.