

Mar 05,
2019

Working

E. coli K12 DNA Extraction

Forked from [E. coli K12 DNA Extraction](#)

[Kenneth Schackart](#)¹, [Kattika Kaarj](#)¹

¹University of Arizona

[dx.doi.org/10.17504/protocols.io.yujfwun](https://doi.org/10.17504/protocols.io.yujfwun)

481b Laboratory

 [Kenneth Schackart](#) 

ABSTRACT

How to extract DNA from *E. coli* K12 using Wizard® Genomic DNA Purification Kit by Promega®.

I do not claim any credit for the development of this protocol. It has been adapted from the protocol detailed in:

 [Wizard Genomic DNA Purification.pdf](#)

PROTOCOL STATUS

Working

We use this protocol in our group and it is working

GUIDELINES

Lab coat and gloves must be worn at all times.

MATERIALS

NAME

Wizard(R) Genomic DNA Purification Kit

CATALOG

A1620

VENDOR

[Promega](#)

MATERIALS TEXT

Reagents in kit:

- Nuclei Lysis Solution
- RNase Solution
- Protein Precipitation Solution
- DNA Rehydration Solution

Additional Reagents:

- Isopropanol
- 70% Ethanol

Additional Materials:

- 2 mL centrifuge tubes
- Pipettes and tips
- Heating block
- Centrifuge

BEFORE STARTING



Spray work area with 70% EtOH solution.

Culture bacteria

- 1 Culture *E. coli* K12 in BHI broth overnight.

 2 mg lyophilized *E. coli* K12 in  10 ml BHI broth.

Pellet the cells



- 2 Add  1 ml cell suspension to 2 mL microcentrifuge tube.
- 3 Label centrifuge tube with your group number.
- 4 Centrifuge at $13,000\text{--}16,000 \times g$ for  00:02:00 .
- 5 Remove supernatant.

Lyse nuclei

- 6 Add  600 μl of Nuclei Lysis Solution.



Nuclei Lysis Solution is marked "NL"




- 7 Gently pipet until the cells are resuspended.
- 8 Incubate at  80 °C on heating block for  00:05:00 to lyse the cells.
- 9 Cool to room temperature.

Degrade RNA

- 10 Add  3 μl RNase Solution to the cell lysate.



RNase solution is marked "R"

- 11 Invert 2-5 times to mix.
- 12 Incubate at  37 °C for  00:15:00 to  01:00:00 .
- 13 Cool to room temperature.


Precipitate proteins

- 14 Add  200 μ l of Protein Precipitation Solution to the RNase-treated cell lysate.




Protein Precipitation solution is marked "P"

- 15 Vortex vigorously at high speed for  00:00:20 .

- 16 Incubate on ice for  00:05:00

- 17 Centrifuge at 13,000-16,000 $\times g$ for  00:03:00 .

Harvest DNA

- 18 Transfer the supernatant containing the DNA to a clean 1.5 mL microcentrifuge tube containing  600 μ l isopropanol.



Some supernatant may remain in the original tube containing the protein pellet. Leave this residual to avoid contaminating the DNA solution with the precipitated protein.



Isopropanol is marked "IPA"


- 19 Label centrifuge tube with your group number.

- 20 Gently mix by inversion until the thread-like strands of DNA form a visible mass.

Wash and dry DNA

- 21 Centrifuge at 13,000-16,000 $\times g$ for  00:02:00 .

- 22 Carefully pour off the supernatant and drain the tube on clean absorbent paper.

- 23 Add  600 μ l of 70% ethanol and gently invert the tube several times to wash the DNA pellet.



70% ethanol solution is marked "Et"

24 Centrifuge at $13,000\text{--}16,000 \times g$ for  00:02:00 .

25 Carefully pour off the ethanol.

26 Drain the tube on clean absorbent paper and allow to air-dry for 10-15 minutes.

Rehydrate DNA

27 Add  100 μl of DNA rehydration solution to the tube.



DNA Rehydration solution is marked "DR"

28 Rehydrate by incubating the solution overnight at room temperature or  4 °C .

29 Store DNA at  2 °C to  8 °C .



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited