

# Total Chlorophyll a Measurements by Spectrophotometer

Dr. Steven Wilhelm

## Abstract

Please contact Dr. Steven Wilhelm (wilhelm@utk.edu) for additional information regarding this protocol.

Adapted from Wetzel and Likens 2000. Limnological Analyses, Springer NY

**Citation:** Dr. Steven Wilhelm Total Chlorophyll a Measurements by Spectrophotometer. **protocols.io**  
dx.doi.org/10.17504/protocols.io.ijpccmn

**Published:** 21 Jun 2017

## Protocol

### Sample Preparation

#### Step 1.

Obtain a sample to be examined. This can be a culture or a field sample.

#### 📌 NOTES

**Alyssa Alsante** 21 Jun 2017

For chlorophyll extractions, it is important to process at least duplicate samples, but triplicate samples are even better.

#### Step 2.

Collect materials on 0.2  $\mu$ M polycarbonate membrane filters with gentle filtration

#### 📌 NOTES

**Alyssa Alsante** 21 Jun 2017

Record the volume that is being filtered in L. The volume to be filtered is dependent on the sample in question. Typically, it will range from 25 mLs to 100 mLs. The more productive (eutrophic) the system, the less you will filter. You should avoid adding too much material as this will clog the filter and it will no longer select based on size exclusion.

#### Step 3.

Place the filter containing the sample into a 15 mL Falcon tube

#### Step 4.

Add 90% acetone to the tube

 **AMOUNT**

5 ml Additional info:

#### Step 5.

Store sample and acetone at -20°C overnight to allow for extraction of Chl a

 **NOTES**

**Alyssa Alsante** 21 Jun 2017

If -20 degree C is not available, 4 degree C will work. The main point is that the sample should be stored in the dark.

Determine Chl a Concentration in Extracts

#### Step 6.

Read the standards in a spectrophotometer at the following absorbance: 665 nm, 645 nm, 730 nm

#### Step 7.

Equation: mg/L Chl a in the extract =  $11.75 \cdot (A_{665} - A_{730}) - 1.31 \cdot (A_{645} - A_{730})$

 **NOTES**

**Alyssa Alsante** 21 Jun 2017

This equation gives results in mg/L, not ug/L.

#### Step 8.

To convert back to the original sample, use the following equation:

Chlorophyll in sample (mg/L) \* acetone volume = total chlorophyll extracted

Total chlorophyll extracted (mg) / volume filtered (L) = concentration of chlorophyll in the original sample (mg/L)