

E. coli Heat Shock Transformation

Snehadri Sinha

Abstract

Citation: Snehadri Sinha E. coli Heat Shock Transformation. **protocols.io**

[dx.doi.org/10.17504/protocols.io.gzcbx2w](https://doi.org/10.17504/protocols.io.gzcbx2w)

Published: 17 Jan 2017

Protocol

Step 1.

Thaw competent cells on ice (20-30min).

Step 2.

Combine 1-5µl of DNA (10pg - 100ng) into 50µL of competent cells in a microcentrifuge tube. Flick the bottom of the tube with your finger a couple of times to mix.

Step 3.

Incubate cell/DNA mixture on ice for 20-30min.

Step 4.

Place transformation tube into a 42°C water bath for 42 seconds (30-60 sec).

Step 5.

Return tube to ice for 2 min.

Step 6.

Add 500µl LB medium, grow in 37°C shaking incubator for 1 h.

Step 7.

Plate transformation on LB agar plates containing the appropriate antibiotic. Often e.g. 50 uL on one plate and 200 uL on another gives a good chance of single colonies. The rest of the transformation can be left on the benchtop overnight and plated the next day if needed.

Step 8.

Incubate plates at 37°C overnight.