



## Preparing primary T cells for fluorescence microscopy

Bulent Arman Aksoy<sup>1</sup>, PINAR AKSOY<sup>1</sup>, Jeff Hammerbacher<sup>1</sup>

<sup>1</sup>Medical University of South Carolina

dx.doi.org/10.17504/protocols.io.vede3a6



TAGS

# microscopy

dapi

Show tags

PROTOCOL STATUS

# Working

We use this protocol in our group and it is working

#### MATERIALS

NAME Y PBS	CATALOG # V	<b>VENDOR</b> ✓ Invitrogen - Thermo Fisher
$BD\ Cytofix/Cytoperm^{\rm TM}\ Fix at ion/Permeabilization\ Solution\ Kit$	554714	BD Biosciences
Pierce™ 16% Formaldehyde (w/v) Methanol-free	28906	Thermo Fisher Scientific
Alexa Fluor™ 488 Phalloidin	A12379	
ProLong™ Glass Antifade Mountant with NucBlue™ Stain	P36983	Thermo Fisher Scientific
Cytology Funnels for Shandon CytoSpin™ Centrifuge Biomedical Polymers	80094-254	VWR international Ltd

# Cell Culture

1 Culture primary T cells using standard tissue culture techniques

## Fixation, permeabilization, and staining

- 2 Collect 200K to 1 million cells and sping them at 300g for  $\boxed{\circlearrowleft 00{:}05{:}00}$
- 3 Discard the supernatant, resuspend the cells in 1 mL PBS, spin them at 300g for © 00:05:00 , and discard the supernatant.
- 4 Fix the cells:

If there are no cell clumps, simply resuspend cells in under the control of Cytofix solution.

© 00:30:00 on a nutating mixer, e.g. Incubate the tubes at 8 4 °C **EQUIPMENT** Fisherbrand™ Nutrating Mixers - Variable Speed Mixer/Shaker Fisherbrand 88-861-043 👄 Spin the cells at 300g for 00:05:00 , discard the supernatant, resuspend them in of 1X CytoPerm Wash Buffer. Optionally, add your stains, dyes (e.g. phalloidin stain) into the permeabilization buffer. § 25 °C (Room temperature) for Cover the tubes with aliminum foil and incubate at **© 00:30:00** or for the specific duration required for the stain on a nutating mixer, e.g. **DEQUIPMENT** Fisherbrand<sup>TM</sup> Nutrating Mixers - Variable Speed Mixer/Shaker Fisherbrand 88-861-043 Spin the cells at 300g for () 00:05:00 , discard the supernatant, resuspend them in 10 Spin the cells at 300g for \( \int \text{00:05:00} \) , discard the supernatant, resuspend them in **200** μl Slide preparation 11 200 µl of the stained cell suspension onto a glass slide or sample cover slip at 300-500g for © 00:05:00 Mount the sample in between a glass slide and cover slip using one drop (at least 10 uL) of ProLong Anti-fade reagent and keep the slide 12 within a dark container for at least \ \( \bigcirc 12:00:00 \)

This is an open access protocol distributed under the terms of the Creative Commons Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited