TRIS acetate phosphate (TAP) medium

Steven Burgess

Abstract

TAP medium is the most widely used for growth of the green alga *Chlamydomonas reinhardtii*. Note the presence of acetate effects the photosynthetic ability of the alga and minimal medium should be used instead if quantifying effects of mutations or inhibitors on photosynthesis.

Source:

from Gorman, D.S., and R.P. Levine (1965) Proc. Natl. Acad. Sci. USA 54, 1665-1669.

http://www.chlamycollection.org/TAP.html

Citation: Steven Burgess TRIS acetate phosphate (TAP) medium. protocols.io

dx.doi.org/10.17504/protocols.io.e95bh86

Published: 23 Jul 2016

Guidelines

From Gorman, D.S., and R.P. Levine (1965) Proc. Natl. Acad. Sci. USA 54, 1665-1669.

This is probably the most widely-used medium at present for experimental work.

Protocol

Step 1.

TRIS base



2 g Additional info:

Step 2.

TAP Salts

■ AMOUNT 25 ml Additional info: **№** PROTOCOL **CONTACT: Steven Burgess**

. TRIS Acetate Phosphate Salts (x40)

Step 2.1.

Ammonium chloride

 NH_4CI

■ AMOUNT

15 g Additional info:

Step 2.2.

Magnesium sulphate

 $MgSO_4$. $7H_2O$

AMOUNT

4 g Additional info:

Step 2.3.

Calcium chloride

CaCl₂ . 2H₂O

■ AMOUNT

2 g Additional info:

Step 2.4.

Distilled H₂O

■ AMOUNT

1 L Additional info:

Step 3.

Phosphate Solution

₽ PROTOCOL

. Phosphate solution

CONTACT: Steven Burgess

Step 3.1.

Potassium hydrogen phosphate

K₂HPO₄ **■** AMOUNT 29 g Additional info: Step 3.2. Potassium dihydrogen phosphate KH₂PO₄ **■** AMOUNT 14 g Additional info: Step 4. Glacial acetic acid **■** AMOUNT 1 ml Additional info: Step 5. **Hutner's Trace Elements ■** AMOUNT 1 ml Additional info: Step 6. Distilled dH₂O NOTES Steven Burgess 21 Jul 2016 Add up to 1L Step 7. Autoclave

NOTES

Steven Burgess 21 Jul 2016

For solid medium add 15g of agar per L.