

10x Tris-Glycine PAGE Running Buffer

Christopher Bartley

Abstract

Western blot running buffer.

Citation: Christopher Bartley 10x Tris-Glycine PAGE Running Buffer. **protocols.io**

[dx.doi.org/10.17504/protocols.io.cq6vzd](https://doi.org/10.17504/protocols.io.cq6vzd)

Published: 31 Jan 2015

Protocol

Step 1.

Fill 1L pyrex bottle with 700mL dH₂O

Step 2.

Add 30.2g Tris base

Step 3.

Add 144.2g glycine

Step 4.

pH solution to 8.80 after dissolution of tris and glycine

Step 5.

Add 10g SDS (1% final)

Step 6.

Fill to 1L with dH₂O