

Chu's-10 Media

Dr. Steven Wilhelm

Abstract

Please contact Dr. Steven Wilhelm (wilhelm@utk.edu) for additional information regarding this protocol.

Citation: Dr. Steven Wilhelm Chu's-10 Media. **protocols.io**

dx.doi.org/10.17504/protocols.io.h8mb9u6

Published: 01 Jun 2017

Protocol

Step 1.

Add 1 L Milli-Q H₂O to a clean media bottle

Step 2.

Add 0.232 g Ca(NO₃)₂*4H₂O



REAGENTS



calcium nitrate by Contributed by users

Step 3.

Add 0.01 g K₂HPO₄



REAGENTS



Potassium phosphate (dibasic) [View](#) by [P212121](#)

Step 4.

Add 0.025 g MgSO₄*7H₂O



REAGENTS



Magnesium Sulfate [View](#) by [P212121](#)

Step 5.

Add 0.02 g Na₂CO₃



REAGENTS



Sodium carbonate [View](#) by [P212121](#)

Step 6.

Add 0.044 g $\text{Na}_2\text{SiO}_3 \cdot 5\text{H}_2\text{O}$



REAGENTS



Sodium metasilicate by Contributed by users

Step 7.

Add 3.5 mg Ferric citrate

Step 8.

Add 3.5 mg Citric acid



REAGENTS



Citric Acid [View](#) by [P212121](#)

Step 9.

Add 1 mL of Metal solution



PROTOCOL

. [Chu's-10 Metal Solution](#)

CONTACT: [Steven Wilhelm](#)

Step 9.1.

Add 1 L Milli-Q H_2O to a clean media bottle

Step 9.2.

Add 2.4 g H_3BO_3



REAGENTS

Boric acid BP1681 by [Fisher Scientific](#)

Step 9.3.

Add 1.4 g $\text{MnCl}_2 \cdot 4\text{H}_2\text{O}$



REAGENTS

Manganese chloride 7773-01-5 by [Fisher Scientific](#)

Step 9.4.

Add 0.4 g ZnCl_2



REAGENTS

✓ Zinc dichloride by Contributed by users

Step 9.5.

Add 0.02 g $\text{CoCl}_2 \cdot 6\text{H}_2\text{O}$



REAGENTS

✓ Cobaltous chloride hexahydrate by Contributed by users

Step 9.6.

Add 0.1 mg $\text{CuCl}_2 \cdot 2\text{H}_2\text{O}$



REAGENTS

Copper (II) chloride dihydrate 10125-13-0 by [Fisher Scientific](#)

Step 10.

If making agar, add 15 g to media

Step 11.

Autoclave at 121°C for 20 min