

# Adult human lung cell dissociation (on ice)

Andrew Potter

## Abstract

This protocol was used to generate a single cell suspension from adult human lung tissue. The procedure is carried out on ice, reducing artifact gene expression changes.

**Citation:** Andrew Potter Adult human lung cell dissociation (on ice). **protocols.io**

[dx.doi.org/10.17504/protocols.io.rcad2se](https://dx.doi.org/10.17504/protocols.io.rcad2se)

**Published:** 27 Jun 2018

## Guidelines

### Enzyme Mixes

#### **Coll. A/Elastase/Dispase Enzyme Mix (make two tubes - each 1 mL)**

60 µL Collagenase A 100 mg/mL – 6 mg/mL final (Sigma, 10103578001)

100 µL elastase 43 u/mL - 4.3 u/mL final (Worthington, LS002292)

100 µL Dispase 90 u/mL – 9 u/mL final (Worthington, LS02100)

5 µL 1 M CaCl<sub>2</sub> – 5 mM final

5 µL DNase (125 U/mL)

730 µL DPBS (no Ca, no Mg)

**+13 mg tissue per 1 mL enzyme mix**

#### BEFORE STARTING

-Prepare enzyme mixes and leave on ice.

-Cool centrifuges to 4 °C.

## Protocol

### **Step 1.**

Transport tissue in ice-cold PBS.

### **Step 2.**

Mince tissue on petri dish on ice using razor blade for 2 min into 1-mm<sup>3</sup> pieces.

### **Step 3.**

Weigh out 13 mg tissue. Using a sterile razor blade or forceps place 13 mg tissue in 1 mL enzyme mix in 1.5 mL eppendorf tube, incubating on ice.

#### AMOUNT

13 mg Additional info:

minced human lung tissue

---

**Step 4.**

Incubate on ice. Triturate 10x using 1 mL pipet set to 700  $\mu$ L every 3 min (w/tip cut). Shake 3-5X to re-suspend every 2 min.

---

**Step 5.**

After 45 minutes of incubation let settle on ice 1 min.

---

**Step 6.**

Remove 80% of the supernatant (consisting of released cells), leaving undigested tissue chunks on the bottom of the tube.

---

**Step 7.**

Add released cells to sterile 30  $\mu$ M filter on 50 mL conical. Rinse filter w/15 mL ice-cold PBS/BSA 0.04%.

 **AMOUNT**

15 ml Additional info: ice-cold PBS/BSA 0.04%

---

**Step 8.**

Divide flow-through into two 15 mL conicals. Bring the volume for each to 14 mL with ice-cold PBS/BSA 0.04%.

 **AMOUNT**

14 ml Additional info: ice-cold PBS/BSA 0.04%

---

**Step 9.**

Spin the two 15 mL conicals with released cells 650 g for 5 min at 4 °C.

 **TEMPERATURE**

4 °C Additional info:

---

**Step 10.**

Remove supernatant for the 15 mL conicals with released cells. Re-suspend the pellets in 14 mL ice-cold PBS/BSA 0.04% for each tube.

 **AMOUNT**

14 ml Additional info: ice-cold PBS/BSA 0.04%

---

**Step 11.**

Add additional 1 mL enzyme mix to residual tissue chunks.

 **AMOUNT**

1 ml Additional info: enzyme mix

---

**Step 12.**

Continue incubating on ice for 35 additional minutes (1 hr. 20 min. total). Triturate 10x using 1 mL pipet set to 700  $\mu$ L every 5 min (w/tip cut). Shake 3-5X to re-suspend every 3 min.

---

**Step 13.**

After 1 hr. 20 min total incubation time triturate digest mix 10X and add digest mix to a new sterile 30  $\mu$ M filter on 50 mL conical.

---

**Step 14.**

Rinse filter w/10 mL ice-cold PBS/BSA 0.04%. Transfer flow-through to two 15 mL conicals.

 [AMOUNT](#)

10 ml Additional info: ice-cold PBS/BSA 0.04%

---

**Step 15.**

Bring the volume for each conical to 14 mL w/ice-cold PBS/BSA 0.04%.

---

**Step 16.**

Spin all four 15 mL conical tubes including the two from the previous step, 650 g for 5 min. at 4 °C.

 [TEMPERATURE](#)

4 °C Additional info:

---

**Step 17.**

Remove supernatant for all tubes. Re-suspend combined volume in 5 mL RBC lysis buffer. Pipet 20x to mix. Let incubate on ice 5 min.

 [AMOUNT](#)

5 ml Additional info: RBC lysis buffer

---

**Step 18.**

Add 10 mL ice-cold PBS/BSA 0.04% to 5 mL RBC lysis buffer. Triturate and apply to sterile 30  $\mu$ M filter on 50 mL conical.

 [AMOUNT](#)

10 ml Additional info: ice-cold PBS/BSA 0.04%

---

**Step 19.**

Transfer flow-through to two 15 mL conicals. Bring volume for each to 14 mL with ice-cold PBS/BSA 0.04%

 [AMOUNT](#)

15 ml Additional info: ice-cold PBS/BSA 0.04%

---

**Step 20.**

Spin 650 g for 5 min at 4 °C. Remove supernatant.

---

**Step 21.**

Re-suspend cells in 250 µL ice-cold PBS/BSA 0.04%. Analyze viability and cell yield using a hemocytometer with trypan blue.

 AMOUNT

250 µl Additional info: ice-cold PBS/BSA 0.04%

---

**Step 22.**

Adjust cell concentration to 1000 cells / µL for 10x chromium or 100 cells/µL for DropSeq.

---