



## Field sampling of root-associated microbes for DNA/RNA extraction

Version 2

Roey Angel<sup>1</sup>

<sup>1</sup>Soil and Water Research Infrastructure

[dx.doi.org/10.17504/protocols.io.vd5e286](https://dx.doi.org/10.17504/protocols.io.vd5e286)

SoWa Anaerobic and Molecular Microbiology



Roey Angel

Soil and Water Research Inrast...



### ABSTRACT

This protocol describes a procedure for sampling plant roots in the field for future DNA and RNA extraction for microbiome analysis. The protocol is deliberately designed to be simple and requires no electronic equipment. Root samples are preserved in LifeGuard Soil Preservation Solution for protecting against nucleic acid degradation.

### TAGS

dna

plant

Show tags

### PROTOCOL STATUS

**Working**

We use this protocol in our group and it is working

### MATERIALS

NAME ▾	CATALOG # ▾	VENDOR ▾
Micro-spatula set	AT16.1	Carl Roth
LifeGuard Soil Preservation Solution	12868-100	Qiagen
Scissors	HCT7.1	Carl Roth
Technical-grade ethanol (70%)	T913.1	Carl Roth
Paper towels	Y03.1	Carl Roth
Microcentrifuge tubes 2 ml	CK06.1	Carl Roth
Garden trowel	View	Amazon
Disposable pasteur pipettes	EA61.1	Carl Roth
Tweezers set	PX40.1	Carl Roth
Cooling box	AA46.1	Carl Roth
Cooling packs	E447.1	Carl Roth

### STEPS MATERIALS

NAME ▾	CATALOG # ▾	VENDOR ▾
Garden trowel	View	Amazon
Micro-spatula set	AT16.1	Carl Roth
Scissors	HCT7.1	Carl Roth


NAME ▾	CATALOG # ▾	VENDOR ▾
Paper towels	<a href="#">Y03.1</a>	<a href="#">Carl Roth</a>
Scissors	<a href="#">HCT7.1</a>	<a href="#">Carl Roth</a>
Paper towels	<a href="#">Y03.1</a>	<a href="#">Carl Roth</a>
Tweezers set	<a href="#">PX40.1</a>	<a href="#">Carl Roth</a>
Microcentrifuge tubes 2 ml	<a href="#">CK06.1</a>	<a href="#">Carl Roth</a>
Micro-spatula set	<a href="#">AT16.1</a>	<a href="#">Carl Roth</a>
LifeGuard Soil Preservation Solution	<a href="#">12868-100</a>	<a href="#">Qiagen</a>
Disposable pasteur pipettes	<a href="#">EA61.1</a>	<a href="#">Carl Roth</a>
Cooling box	<a href="#">AA46.1</a>	<a href="#">Carl Roth</a>
Cooling packs	<a href="#">E447.1</a>	<a href="#">Carl Roth</a>

#### BEFORE STARTING

- Clean spatulas using 70% ethanol


**1 Sample a triplicate of plant individuals spaced out a few metres apart from each other** Make sure you are really sampling individual plants and not offshoots of the same plant.

**2 Using a garden trowel, carefully dig out the plants while keeping the root system intact** (as much as possible, of course).




Garden trowel  
by [Amazon](#)  
[View](#)

**3 While holding the plant by the shoot, shake the root system hard enough so that all loose soil is removed from it.** Take care to damage the plant as little as possible You can use a spatula to remove large soil aggregates that are attached to the roots.




Micro-spatula set  
by [Carl Roth](#)  
Catalog #: [AT16.1](#)

**4 From the remaining root system (plus soil particles plus attached to the roots), trim a 'representative sample' of roots using scissors or scalpel.** It is usually best to trim the roots onto a piece of paper towel.



Scissors  
by [Carl Roth](#)  
Catalog #: [HCT7.1](#)




Paper towels

by [Carl Roth](#)

Catalog #: [Y03.1](#)


**5 Cut the trimmed out roots a little so that they fit into a 2.0 ml tube.**



Scissors

by [Carl Roth](#)

Catalog #: [HCT7.1](#)




Paper towels

by [Carl Roth](#)

Catalog #: [Y03.1](#)


**6 Place about 2-3 g of that cut out sample into a 2.0 ml tube.**



Tweezers set

by [Carl Roth](#)

Catalog #: [PX40.1](#)



Microcentrifuge tubes 2 ml

by [Carl Roth](#)


Catalog #: [CK06.1](#)



**NOTE**

The root tissue should make up at least half or more of the mass, while the remaining attached soil should make up the rest

**7 Press the sample a little into the bottom of the tube to decrease its volume.**



Micro-spatula set


by [Carl Roth](#)

Catalog #: [AT16.1](#)


**NOTE**

Make sure the roots do not take up more than 1/2-2/3 of the volume  
It is of course possible to split each sample into several separate tubes, depending on the specific type of roots, and submerge each with LifeGuard solution

**8 Add as much LifeGuard solution so that the sample is submerged in about twice of its volume (about 1.0 – 1.5 ml).**  
Best is to use a disposable Pasteur pipette for dispensing the solution.




LifeGuard Soil Preservation Solution  
by [Qiagen](#)  
Catalog #: [12868-100](#)




Disposable pasteur pipettes  
by [Carl Roth](#)  
Catalog #: [EA61.1](#)

 1.5 ml

- 9 Place the tubes in cooling (around 4 °C) and keep them cooled until you reach the lab. The solution will protect nucleic acids even at room temperature for several days, but cooling is preferred.



Cooling box  
by [Carl Roth](#)  
Catalog #: [AA46.1](#)




Cooling packs  
by [Carl Roth](#)  
Catalog #: [E447.1](#)

 4 °C

- 10 In the lab, store the samples in a freezer (-20 – -80 °C).

 -20 °C or -80 °C

 This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited