

Artificial Seawater Based AMP1 Medium

Bonnie Poulos

Abstract

Citation: Bonnie Poulos Artificial Seawater Based AMP1 Medium. **protocols.io**

dx.doi.org/10.17504/protocols.io.c9mz45

Published: 04 Jan 2016

Guidelines

Turk's Island Salt Mix

Chemical	g/2 L	Final Conc.	Mfc.-PN	Chemical Grade
NaCl	56.22	481 mM	Fisher S271	ACS
MgSO ₄ ·7H ₂ O	13.8	28 mM	Sigma M2773	Mo.Biol. 99%
MgCl ₂ ·6H ₂ O	10.98	27 mM	AlfaAesar	99.99%
CaCl ₂ ·2H ₂ O	2.94	10 mM	Fisher C79	ACS
KCl	1.34	9 mM	Acros	99+%
MQ·H ₂ O	QS to 2000 ml	na	in house	na

Macronutrients

Macronutrient	g/100 ml	Final Conc.	Mfc.-PN	Chemical Grade
0.1 M NaH ₂ PO ₄ ·H ₂ O	1.38	50 µM	Mallinckrodt 7892-04	ACS
0.8 M (NH ₄) ₂ SO ₄	10.57	400 µM	VWR-BDH 0216	ACS 99%

Buffers

Buffer	Stock	Final Conc.	Mfc.-PN	Chemical Grade
0.6 M NaHCO ₃	0.504 g/10 ml	6 mM	Sigma S6014	ACS 99.7-100.3%
1 M HEPES	11.9 g/50 ml	1 mM	Fisher BP310	Mol. Biol.


Trace Metal Mix Working Stock

Trace Metal	Primary Stock	Final Conc.	Mfc.-PN	Chemical Grade
Na ₂ EDTA·2H ₂ O	0.2175 g/50 ml	0.1170 µM	Sigma E4884	ACS 99-101%
FeCl ₃ ·6H ₂ O	0.16 g/50 ml	0.1180 µM	Sigma 44944	ACS 98-102%
ZnSO ₄ ·7H ₂ O	1.15 g/50 ml	0.0008 µM	Sigma 204986	>99.5%
CoCl ₂ ·6H ₂ O	0.595 g/50 ml	0.0005 µM	Sigma 60820	>98%
MnCl ₂ ·4H ₂ O	8.905 g/50 ml	0.0090 µM	M3634	ReagentPlus >99%
Na ₂ MoO ₄ ·2H ₂ O	0.363 g/50 ml	0.0003 µM	M1651	>99.5%
Na ₂ SeO ₃	0.865 g/50 ml	0.0010 µM	S5261	BioReagent ~98%

NiCl₂·6H₂O 1.19 g/50 ml 0.0010 μM 223387 ReagentPlus

Note: Use primary stocks of ZnSO₄, CoCl₂, MnCl₂, Na₂MoO₄, Na₂SeO₃, and NiCl₂ prepared for Pro99 trace metal working stock.

Math Calculations				
Chemical	Molar Mass (g/mol)	Amount in Amp-1	Final Conc	
NaCl	58.44 g/L	56.23 g/L	$\left(\frac{56.23}{2000}\right) \times \left(\frac{58.44}{1000}\right)$	$= 0.481 \text{ M} = 481 \text{ mM}$
MgSO ₄ ·7H ₂ O	246.47 g/L	13.8 g/L	$\left(\frac{13.8}{2000}\right) \times \left(\frac{246.47}{1000}\right)$	$= 0.02799 \text{ M} = 28 \text{ mM}$
MgCl ₂ ·6H ₂ O	203.30 g/L	10.98 g/L	" "	$= 0.02700 \text{ M} = 27 \text{ mM}$
CaCl ₂ ·2H ₂ O	147.01 g/L	2.94 g/L	" "	$= 0.00999 \text{ M} = 10 \text{ mM}$
KCl	74.55 g/L	1.34 g/L	" "	$= 0.00898 \text{ M} = 9 \text{ mM}$
NaH ₂ PO ₄ ·H ₂ O	137.99 g/L	$\left(\frac{0.25 \text{ ml}}{250 \text{ ml}}\right) \times \left(\frac{1.28 \text{ g}}{100 \text{ ml}}\right)$	" "	$= 0.00054 = 0.54 \text{ mM} = 54 \mu\text{M}$
(NH ₄) ₂ SO ₄	132.14 g/L	$\left(\frac{0.25 \text{ ml}}{250 \text{ ml}}\right) \times \left(\frac{10.57 \text{ g}}{100 \text{ ml}}\right)$	" "	$= 0.00041 \text{ M} = 0.41 \text{ mM} = 409 \mu\text{M}$
NaHCO ₃	84.01 g/L	$\left(\frac{5 \text{ ml}}{250 \text{ ml}}\right) \times \left(\frac{0.504 \text{ g}}{10 \text{ ml}}\right)$	" "	$= 0.00599 \text{ M} = 6 \text{ mM}$
HEPES	238.3 g/L	$\left(\frac{0.0105 \text{ ml}}{250 \text{ ml}}\right) \times \left(\frac{11.9 \text{ g}}{50 \text{ ml}}\right)$	" "	$= 0.0009987 = 1 \text{ mM}$
Na ₂ EDTA·2H ₂ O	372.24 g/L	$\frac{0.005 \text{ ml}}{500 \text{ ml}} \times \frac{0.245 \text{ g}}{50 \text{ ml}}$	" "	$= 0.000001 = 0.1 \mu\text{M}$
FeCl ₃ ·6H ₂ O	270.3 g/L	$\frac{0.009 \text{ ml}}{500 \text{ ml}} \times \frac{0.16 \text{ g}}{50 \text{ ml}}$	" "	$= 1.8392 \times 10^{-7} = 0.18 \mu\text{M}$
ZnSO ₄ ·7H ₂ O	287.56 g/L	$\frac{0.005 \text{ ml}}{1 \times 10^{-5} \times 500 \text{ ml}} \times \frac{50 \mu\text{M}}{1.13 \text{ g/50 ml}} \times (0.05 \text{ ml})$	" "	$= 3.9932 \times 10^{-8} =$
CoCl ₂ ·6H ₂ O	237.93 g/L	$\frac{0.005 \text{ ml}}{1 \times 10^{-5} \times 500 \text{ ml}} \times \frac{50 \mu\text{M}}{1.19 \text{ g/50 ml}} \times (0.05 \text{ ml})$	" "	$= 2.5007 \times 10^{-8} = 5 \times 10^{-9} \text{ M}$
<div style="border: 1px solid black; border-radius: 50%; padding: 10px; display: inline-block;">See next page</div>				
$ \begin{aligned} & \frac{2.5493 \times 10^{-4}}{1 \times 10^{-5}} \times \frac{1.19 \times 10^{-4}}{1.19 \times 10^{-4}} \\ & \quad \downarrow \\ & \quad 1.19 \times 10^{-9} \\ & \quad \downarrow \\ & \quad 1.19 \times 10^{-9} \\ & \quad \downarrow \\ & \quad 5 \times 10^{-9} \end{aligned} $				



Materials

Ethylenediaminetetraacetic acid disodium salt dihydrate [E4884](#) by [Sigma Aldrich](#)

Sodium Chloride [S271](#) by [Fisher Scientific](#)

Magnesium sulfate heptahydrate [M2773](#) by [Sigma Aldrich](#)

Calcium Chloride Dihydrate [C79](#) by [Fisher Scientific](#)

Sodium bicarbonate [S6014](#) by [Sigma Aldrich](#)

HEPES [BP310](#) by [Fisher Scientific](#)

Iron(III) chloride hexahydrate [44944](#) by [Sigma Aldrich](#)

Zinc sulfate heptahydrate [204986](#) by [Sigma Aldrich](#)

Cobalt(II) chloride hexahydrate [60820](#) by [Sigma Aldrich](#)

Manganese(II) chloride tetrahydrate [M3634](#) by [Sigma Aldrich](#)

Sodium molybdate dihydrate [M1651](#) by [Sigma Aldrich](#)

Sodium selenite [S5261](#) by [Sigma Aldrich](#)

Nickel(II) chloride hexahydrate [223387](#) by [Sigma Aldrich](#)

Protocol

Turk's Island Salt Mix

Step 1.

Dissolve each salt completely before adding the next one

NOTES

VERVE Team 12 Aug 2015

Refer to table in guidelines for full list of salts.

Turk's Island Salt Mix

Step 2.

Dispense into 500 ml acid-washed polycarbonate bottles

Turk's Island Salt Mix

Step 3.

Autoclave 30 min.

DURATION

00:30:00

Macronutrients

Step 4.

pH NaH_2PO_4 to 7.5 using 1M NaOH

NOTES

VERVE Team 01 Jul 2015

About 8ml for 100ml volume.

VERVE Team 01 Jul 2015

Prepare each one separately.

Macronutrients

Step 5.

Filter sterilize each solution using 0.2 μm syringe filter into new, sterile 50ml centrifuge tubes or acid washed and sterile polycarbonate bottles

Macronutrients

Step 6.

Store at 4°C

NOTES

VERVE Team 01 Jul 2015

Dispense 250 μl of each per 500 ml bottle of Turk's Island Salt Mix when preparing final medium.

Buffers

Step 7.

pH HEPES to 7.5 using 1M NaOH

REAGENTS

HEPES [BP310](#) by [Fisher Scientific](#)

NOTES

VERVE Team 01 Jul 2015

About 9 ml for 50 ml volume.

VERVE Team 01 Jul 2015

Prepare each one separately.

Buffers

Step 8.

Filter sterilize each solution using 0.2 µm syringe filter into new, sterile 15 or 50 ml centrifuge tubes or acid washed and sterile polycarbonate bottles

Buffers

Step 9.

Store at 4°C

📌 NOTES

VERVE Team 01 Jul 2015

The NaHCO₃ should be made monthly.

Buffers

Step 10.

Dispense 5 ml of NaHCO₃ and 0.5 ml HEPES per 500 ml bottle of Turk's Island Salt Mix when preparing final medium

Trace Metal Mix Working Stock

Step 11.

Weigh out 0.2175 g Na₂EDTA·2H₂O using dust free paper

📋 AMOUNT

0 g Additional info:

🧴 REAGENTS

Ethylenediaminetetraacetic acid disodium salt dihydrate [E4884](#) by [Sigma Aldrich](#)

Trace Metal Mix Working Stock

Step 12.

Transfer to acid washed 50 ml volumetric flask filled with 40 ml MQ-water

Trace Metal Mix Working Stock

Step 13.

Dissolve EDTA by inverting flask several times

📌 NOTES

VERVE Team 01 Jul 2015

May have to heat 5 min. at 80°C to dissolve.

Trace Metal Mix Working Stock

Step 14.

Weigh out 0.16 g FeCl₃·6H₂O using dust free paper

📋 AMOUNT

0 µl Additional info:

🧴 REAGENTS

Iron(III) chloride hexahydrate [44944](#) by [Sigma Aldrich](#)

Trace Metal Mix Working Stock

Step 15.

Dissolve iron chloride into same volumetric flask by inverting several times

Trace Metal Mix Working Stock

Step 16.

Individually add and dissolve 50 µl each of the ZnSO₄, CoCl₂, MnCl₂, Na₂MoO₄, NaSeO₃, and NiCl₂

Primary Trace Metal Stocks

Trace Metal Mix Working Stock

Step 17.

Adjust volume to 50 ml mark with MQ-water

Trace Metal Mix Working Stock

Step 18.

Filter through a 0.2 µm syringe filter into sterile, acid washed container in laminar flow hood

Trace Metal Mix Working Stock

Step 19.

Store sterile stock at 4°C

Trace Metal Mix Working Stock

Step 20.

Dispense 5 µl Stock Trace Metal Mix to 500 ml bottle of Turk's Island Salt Mix when preparing final medium