



Feb 27, 2019

Working

iPSC Freezing

In 1 collection

Celeste Karch¹, Rita Martinez¹, Jacob Marsh¹¹Washington University in St Louis

dx.doi.org/10.17504/protocols.io.x8zfrx6

Neurodegeneration Method Development Community

Tech. support email: ndcn-help@chanzuckerberg.com

Celeste Karch

Washington University in St Louis



PROTOCOL STATUS

Working

We use this protocol in our group and it is working

GUIDELINES

This protocol is part of the [Screening Edited iPSC Clones collection](#).

STEPS MATERIALS

NAME

CATALOG

VENDOR

Accutase™ Cell Dissociation Reagent

A1110501

Gibco - Thermo Fischer

SAFETY WARNINGS

Please refer to the SDS (Safety Data Sheet) for information about hazards, and to obtain advice on safety precautions.

1 Aspirate media

2 Gently wash cells with 1x PBS.



Use 2-3 mL per well in 6 well plate.

3 Add Accutase (Gibco A11105-01) directly to the cells and incubate at **37 °C** for 3-4 minutes. **00:03:00**

Individual donor cell lines exhibit variable sensitivity to accutase-mediated dissociation. Thus, monitor cells closely to determine when single cell dissociation is achieved.



For a 6 well plate, add 0.75-1 mL per well.

For a 10cm² dish, add 3 mL.



Accutase™ Cell Dissociation Reagent

by [Gibco - Thermo Fischer](#)

Catalog #: [A1110501](#)

4 Tap dish to aid in dislocation of cells.

5 Add DMEM/F12 directly to cells.



- For a 6 well plate, add 2 ml - 4 ml per well.
- For a 10cm² dish, add 9 ml
- If cells remain attached, use a cell scraper to gently dislodge cells (apply gentle pressure and use 1-2 passes to remove cells)

6 Collect cells in conical tube (15mL/50mL depending on volume).

7 Add 2 ml - 5 ml DMEM/F12 to dish to remove all cells from the dish and add to conical tube.

8 Centrifuge cells at 750 rpm for 00:03:00 at room temperature.

9 Carefully aspirate supernatant.



To avoid aspirating cell pellet, it is OK to leave a small amount of media (0.5-1mL).



10 Resuspend cell pellet with mTesR1 (No Rock Inhibitor).

- Use volume appropriate for freezing.
- Assume 1 mL per cryovial total and add ½ total volume of mTesR1.
- Pipet cells 1-2 times only to preserve cell clumps.



Example: to freeze 10 tubes, you will need 10 mL total and will add 5 mL mTesR1 to cell pellet (and 5 mL of 2x Freezing Media below).

11 Add an equal volume of cold 2x Freezing Media (20% DMSO, FBS). Pipet cells 1 time only to preserve cell clumps.

- 12 Transfer cell suspension to pre-labeled cryovials (1 mL per cryovial).
- Ensure that cryovials are labeled with the following
 - Cell Type
 - Line Name
 - Passage #
 - Date
 - Your Name
- 13 Freeze vials at  **-80 °C** in foam racks for 48-72 hours.  **48:00:00**
- 14 Transfer vials to liquid nitrogen for long-term storage.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited