



Mar 28,  
2019

Working

## Chelex DNA isolation for quick plant genotyping

Magdalena Julkowska<sup>1</sup>

<sup>1</sup>King Abdullah University of Science and Technology

[dx.doi.org/10.17504/protocols.io.pdudi6w](https://doi.org/10.17504/protocols.io.pdudi6w)



Magdalena Julkowska  
King Abdullah University of Science and Technology



### ABSTRACT

Quick and dirty, but very very fast protocol for DNA isolation. It works beautifully if you have plenty of lines to genotype - example T-DNA insertion lines for validation of your favourite gene.

### PROTOCOL STATUS

#### Working

We use this protocol in our group and it is working

### MATERIALS

NAME ▾	CATALOG # ▾	VENDOR ▾
Chelex 100	C7901-100G	Sigma Aldrich

### SAFETY WARNINGS

- 1 Grind tissue in liquid N<sub>2</sub> (1 leaf should be enough)

If you collect the tissue in a tube containing glass beads (1-2 mm diameter) you can put the frozen tissue samples in the tissue grinder. This is by far the most efficient method to grind large quantities of various samples.

- 2 Add 200µl 10% Chelex (in MilliQ)

- 3 Vortex

- 4 Incubate 15min @ 95 degrees in shaker

- 5 Vortex

- 6 Spin 15min @ max RPM

- 7 Use (2µl) supernatant for PCR (optional: transfer supernatant to fresh tube for storage at -20C)

Genotyping PCR (DreamTaq)

8

step 1 94 degrees 4min  
step 2 94 degrees 45s  
step 3 50 degrees 45s (try higher temp for more specificity)  
step 4 72 degrees 1:20min (depending on size insert)  
step 5 repeat 35x step 3  
step 6 72 degrees 10min  
step 7 16 degrees 10min

Run the gel

9

Run 20uL of the PCR product on the 1% agarose gel



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited