

# Transient transfection of unicellular relative to animals Abeoforma whisleri with electroporation using Neon System

Maria Rubio-Brotons

## Abstract

**Citation:** Maria Rubio-Brotons Transient transfection of unicellular relative to animals Abeoforma whisleri with electroporation using Neon System. **protocols.io**

dx.doi.org/10.17504/protocols.io.hmvb466

**Published:** 08 Apr 2017

## Protocol

Centrifuge 5 mL of culture at 1000g 5min

### Step 1.

young culture: 2-3 days passage

Resuspend cells in 5mL Sorbitol 1M

### Step 2.

Count cells: between 1.000 and 10.000 cell

### Step 3.

Mix cells with plasmid DNA (between 5 to 10 ug)

### Step 4.

Neon system only allows 10 uL final volume in the tip-- adjust your cell and DNA dilutions accordingly for a final volume 10 uL

Electroporate cells with plasmid in Neon tip

### Step 5.

Settings:

1200V

25ms

1pulse

Transfer cells into chamber with growth medium for recovery overnight

### Step 6.

Medium of Abeoforma whisleri is Marine Broth, each well: 500 uL medium

Observe transformants on the next day after electroporation

**Step 7.**

Efficiency varies from 1 to 50 transformants per 10.000 cells