



Nov 02,
2019

CRISPR RNP Electroporation Protocol

Aditya Mohan¹

¹Johns Hopkins University

1

Works for me

dx.doi.org/10.17504/protocols.io.8yfhxtn



Aditya Mohan



Culture Cells

- 1 **Subculture cells for a minimum of 2-3 days before electroporation, and visually inspect the cells with a microscope to ensure healthy appearance.**
- 2 Resuspend Oligos in IDTE Buffer at 100 uM
For a 2nmol IDT order, dilute oligos with 20 uL of Tris-HCL (Dharmacon)

Form the RNP Complex

- 3 For each well-undergoing electroporation, dilute the guide RNA and Cas9 enzyme in PBS gently swirling the pipet tip while pipetting:
1.3 uL of STERILE PBS
5.1 uL of **GenCrispr NLS-Cas9-EGFP Nuclease**
3.6 of diluted sgRNA oligos

Total: 10 uL of total volume

- 4 Incubate for 15 minutes at Room Temperature

⌚ 00:15:00

Prepare Nucleofector System

- 5 Turn on the nucleofector and load the X-core and the program of the desired cell type (or closest cell type)

Resuspend the Alt-R Cas9 Electroporation Enhancer

- 6 Resuspend the Alt-R Cas9 Electroporation Enhancer to 100 uM in IDTE.

Prepare the Sorting Plate

- 7 Prepare a 96 well culture plate to receive cells following Nucleofection
(Ideally U Bottom plate and at least 175 uL of culture media) and prewarm to 37 degrees C.

Prepare the Quenching Plate

- 8 Prepare a 24 well plate with at least 1 mL of media and prewarm to 37 degrees C

Collect Cells

- 9 Collect 1M cells per electroporation and MAKE SURE TO WASH THE CELLS WITH PBS.
Centrifuge the cells at 30g for 10 minutes at RT
- 10 Remove as much supernatant as possible without perturbing the pellet.

Electroporation

- 11 Resuspend the cell pellet in 96 uL of the Ingenio® Electroporation media
Add 4 ul of the diluted electroporation enhancer
Add 10 uL of the RNP gRNA mix to the well

Zap Step

- 12 Load solution into the Nucleofector module and make sure there are no air bubbles.
Zap and keep in the zap mixture for 10 seconds before removing

Harvest

- 13 Remove the cells from the cuvette and add the cells to the quenching plate for a minimum of 4-5 hours

FACS Sorting

- 14 Take samples to be sorted.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited