

Human Breast Tissues Dissociation

Quy Nguyen¹

¹University of California, Irvine

dx.doi.org/10.17504/protocols.io.twbepan

Human Cell Atlas Method Development Community

Working

Sep 30, 2018










Quy Nguyen ⚡

PROTOCOL STATUS

Working

We use this protocol in our group and it is working

MATERIALS

	NAME ▾	CATALOG # ▾	VENDOR ▾
	Collagenase Type I powder	17100017	Thermo Fisher Scientific
	DMEM and Hams F-12 50/50	MT10090CV	Fisher Scientific
	Deoxyribonuclease I	D4263-5VL	Sigma Aldrich
	Dulbeccos Phosphate-Buffered Salt Solution 1X	MT21031CV	Fisher Scientific
	0.05% Trypsin	MT25051CI	Fisher Scientific
	Antibiotic-Antimycotic	15240062	Thermo Fisher Scientific
	Fetal Bovine Serum	FB-12	Omega Scientific, Inc.

BEFORE STARTING

Prepare 4mg/mL collagenase solution before start and place on 37C shaker to dissolve.

Prepare materials

1

Collagenase Solution

100mg/mL Collagenase Type I in DMEM/Ham's F12

Media

DMEM/Ham's F12

10% FBS

1% Antibiotic

DMEM/Ham's F12

5% FBS

1% Antibiotic

DNase

1mg/mL DNase in PBS

Initial Tissue Preparation

2

Transfer tissues to 150cm plate.

Wash with large quantities of PBS to remove blood and old storage media.

-Amount of PBS varies depending on volume of tissues washed.

Aspirate off PBS.

Repeat PBS wash twice, for a total of three washes.

Collect small pieces of tissues for FFPE or OCT Blocks.

Mechanical Dissociation of Tissue

3

Use two number 10 scalpel to chop tissue.

Chop tissue until epithelial tissues are about 2-3mm in size.

-Depending on staining tissues sample, removal of excessive adipose tissues will help in the chopping and digesting process.

Enzymatic Dissociation of Tissue

4

Transfer about 20mL in volume of chopped tissue to each 50mL conical tube.

Add 20mL of DMEM with 5% FBS with 4mg/mL Collagenase Type I.

Place on 37°C shaker at about 180-200rpm.

Leave on heated shake for 12-16 hours.

Organoid Collection

5

Centrifuge tubes at 400g x 5min.

Aspirate to remove supernatant.

Wash with 50mL of PBS.

Centrifuge at 400g x 5min to remove PBS.

Single Cell Preparation

6

Add 2mL of 0.05% Trypsin to pellet of organoids.

Place at 37°C for 6 minutes.

Take tube out every 2 minutes to pipette up and down with P1000, for a total of three times.

Add 10mL of DMEM with 10% FBS to tube.

Centrifuge at 400g x 5min.

Aspirate off supernatant.

Resuspend pellet in 1mL of DMEM with 10% FBS

Add 100µL of DNase, incubate for 5min.

Add 10mL of DMEM with 10% FBS.

Centrifuge at 400g x 5min

Remove supernatant

Resuspend in 10mL of DMEM with 10% FBS

Count cells.

Centrifuge at 400g x 5min.

Bring back up to appropriate concentration for FACS staining.

Single Cell FACS Preparation

7

Prepare appropriate tubes for FACS staining.

To all tube add the following antibodies.

-PE-CD49f

-APC-EpCAM

-DAPI-CD45/CD31

--Use antibody concentration suggested by manufacture.

Stain for 20min at 4°C

Take out and wash cells with DMEM.

Pass through FACS tube filter.

Before sort begins add SytoxGreen to tube.

Sort for all live cells, anything that is SytoxGreen negative.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited