



May 17,
2019

Working

Yale - Total Protein [↗](#)

John Stack¹, Gary Cline¹

¹Yale University

[dx.doi.org/10.17504/protocols.io.y4afyse](https://doi.org/10.17504/protocols.io.y4afyse)

Mouse Metabolic Phenotyping Centers
Tech. support email: info@mmpc.org

Lili Liang

ABSTRACT

Summary:

Procedure to determine the total protein concentration in blood, serum, and plasma. Total Protein is determined by the biuret reaction.

EXTERNAL LINK

<https://mmpc.org/shared/document.aspx?id=209&docType=Protocol>

MATERIALS

NAME	CATALOG #	VENDOR
Total Protein Reagent	R84086	Prolabs(cliniqa)
Multi Analyte Calibrator	R60010	Prolabs(cliniqa)
Assayed Control Serum 1	R83082	Prolabs(cliniqa)
Assayed Control Serum 2	R83083	Prolabs(cliniqa)

MATERIALS TEXT

Reagent Preparation:

Total Protein Reagent: As supplied by vendor

Multi Analyte Calibrator: Add the appropriate amount of water (6.5mL) to the chemical control bottle.
Invert to mix, allowing 15 minutes for the reagent to settle.

Assayed Control Serum 1: Add the appropriate amount of water (6.5mL) to the chemical control bottle.
Invert to mix, allowing 15 minutes for the reagent to settle.

Assayed Control Serum 2: Add the appropriate amount of water (6.5mL) to the chemical control bottle.
Invert to mix, allowing 15 minutes for the reagent to settle.

BEFORE STARTING

Analysis by automated system Cobas Mira Plus

1 Calibrate Cobas for Total Protein analysis by running a multi-analyte calibrator and two control serum.

2 Sample handling as performed by Cobas Mira Plus

- a) Pipette 5 μ L of sample into cuvette.
- b) Add 150 μ L of Total Protein Reagent.
- c) Incubate at 37°C for 10 minutes
- d) Absorbance is measured at 550nm



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited