

ViparSpectra LED Aquarium Light Setup

Jack Koch, Virginia Weis

Abstract

Citation: Jack Koch, Virginia Weis ViparSpectra LED Aquarium Light Setup. **protocols.io**

dx.doi.org/10.17504/protocols.io.j5zcq76

Published: 02 Oct 2017

Guidelines

Channel 1: Blue light

Channel 2: White/green/red light

Protocol

Turn ON/OFF

Step 1.

With remote: Press power button

With on-light controls: Press in knob and hold for a few seconds. The word "on" or "off" will appear on screen. Release the knob.

Set time

Step 2.

With remote: Press the Clock button. Press up or down to set hour, then press OK. Press up or down to set minute. Wait until H-SE appears on screen.

With on-light controls: Press knob in several times until H-SE appears on screen. Turn knob until Cloc appears on screen. Press knob in to select hour slot. Turn knob to set hour. Press knob in to select minute slot. Turn knob to set minute. Wait until H-SE appears on screen.

Set light intensity

Step 3.

With remote: Press the menu button until H-SE appears. Press OK until Ch1 appears. Press up or down

until desired light intensity is achieved (0-100%). Press OK until Ch2 appears. Press up or down until desired light intensity is achieved (0-100%). Press OK until H-SE appears. Wait until time is shown on screen.

With on-light controls: Press knob in several times until Ch1 appears on screen. Turn knob until desired light intensity is achieved (0-100%). Press knob in several times until Ch2 appears on screen. Turn knob until desired light intensity is achieved (0-100%). Press knob in several times until H-SE appears. Wait until time is shown on screen.

Set channel on/off times

Step 4.

With remote: Press C1 On button. Press up or down to select hour slot. Set hour time. Press OK to select minute slot. Set minute time. Wait until H-SE appears on screen. Repeat for C1 Off, C2 On, C2 Off.

With on-light controls: Press knob in several times until H-SE appears. Turn knob until C1 On appears. Press knob in to select hour slot. Turn knob to set hour. Press knob in to select minute slot. Turn knob to set minute. Wait until H-SE appears on screen. Repeat for C1 Off, C2 On, C2 Off.