

# Making spike-in transcripts for mRNA normalization Version 2

**Darach Miller**

## Abstract

This is adapted from Benjy Neymotin's implementation of the SP6 in vitro transcription protocol for the purposes of doing 4tU labeled transcripts.

**Citation:** Darach Miller Making spike-in transcripts for mRNA normalization. **protocols.io**

dx.doi.org/10.17504/protocols.io.fgnbjve

**Published:** 06 Sep 2016

## Guidelines

Currently, the polyA pSP64 plasmids in Gresham lab are DGP 104,105,106,229. These are the 700,900,1200, and the c elegans spike ins.

## Protocol

### Preparing templates

#### Step 1.

Grow DH5alpha ecoli cells overnight with pSP64 poly(A) plasmid containing inserted spike-in sequence.

#### Step 2.

Miniprep plasmids, quantify. You'll need a lot of plasmid.

#### Step 3.

Set up the digestion reaction in 50ul 1x EcoRI compatible buffer with 3ug pSP64-NSERT-poly(A) and 5ul EcoRI. Incubate 2 hours each, 37C.

Save some plasmid for the gel below.

## NOTES

**Darach Miller** 29 Aug 2016

Benjy recommends you prepare at least two reactions per spike-in.

#### Step 4.

Pour a 1% agarose gel and run 100ng of digest and undigested plasmid, to check for cutting.

#### Step 5.

Assuming you observe complete digestion, clean up reactions using PCR clean up kit, pooling repeats of the same plasmid digest at this stage. Quantify using nanodrop to have an idea for concentration.

#### 📌 NOTES

**Darach Miller** 29 Aug 2016

Concentration should be in the range of 100-200 ng/ul, when resuspended in 30ul.

#### in-vitro txn with SP6

#### Step 6.

Set up SP6 reaction master mix, from the kit reagents:

	1x	5x
Txn optimized 5x buffer	4ul	20ul
DTT (100mM)	2ul	10ul
Recombinant RNasin	0.75ul	3.75ul
10mM rATP	1ul	5ul
10mM rUTP	1ul	5ul
10mM rCTP	1ul	5ul
10mM rGTP	1ul	5ul
4tUTP (10mM)	2ul	10ul
SP6 RNA Polymerase	1ul	5ul

#### 📌 NOTES

**Darach Miller** 06 Sep 2016

Do this on bench, not on ice (the mix). The manual claims cold will encourage precipitation of your DNA with aid of the spermidine in the buffer.

#### Step 7.

Make new tubes of 6.25ul of linearized plasmid DNA. For a kit positive control, dilute 1ul of the provided standard with 5.25ul of water.

Concentrations should be in the range of 100-200 ng/ul.

### **Step 8.**

Add 13.75ul of SP6 reaction master mix to each tube of linearized DNA. Incubate 1-2hr in 30/37C water bath.

#### **NOTES**

**Darach Miller** 29 Aug 2016

Darach has seen that 2 hours at 37C works for thiolated spikeins, but with multiple bands for the celegans spike-in. 2 hours at 30C works well for non-thiolated spikeins.

The manual says use 30C if you get multiple txn products, 37C or 42C if you need more product.

### **Step 9.**

Remove 2ul into PCR tubes for later gel.

### **Step 10.**

To the rest, add 1ul of RQ1 RNase-Free DNase. Incubate 15min at 37C.

### **Step 11.**

To each reaction, add 40ul Ampure XP beads and mix well with pipette. Let sit RT for 5min.

### **Step 12.**

Collected beads to side in a magnetic rack. Aspirated supernatant, all.

### **Step 13.**

Add 500ul 80% etOH onto beads. Let sit 30 seconds, then aspirate off.

### **Step 14.**

Add 500ul 80% etOH onto beads. Let sit 30 seconds, then aspirate off. Make sure to get everything out of the bottom, small tip may help.

### Step 15.

Let beads dry in the rack with a kimwipe over the open tubes, 10min RT.

### Step 16.


Resuspend beads in 20ul hyclone H<sub>2</sub>O.

### Step 17.

Run 3ul on 1% agarose TAE gel for about 20min at 100V with the NEB 1kb ladder. Visualize with your favorite dye (we use sybrsafe now).

#### EXPECTED RESULTS

1kb ladder	Positive control	Positive control, DNAsed	DGP104 reaction	DGP104 reaction, DNAsed	DGP105 reaction	DGP105 reaction, DNAsed	DGP106 reaction	DGP106 reaction, DNAsed	DGP229 reaction	DGP229 reaction, DNAsed	1kb ladder
------------	------------------	--------------------------	-----------------	-------------------------	-----------------	-------------------------	-----------------	-------------------------	-----------------	-------------------------	------------



( DNAsed means also cleaned up with beads )

### Step 18.

Quantify on qubit, dilute. 0.1ng/ul is a good working mix.

Benjy suggests using 4ng spikein (each) : 100ug total RNA. Darach's had success with 1ng for 5e7 cells (for qPCR normalization).