



BG11 medium (working group Wilde)

Christin Köbler¹, Annegret Wilde¹

¹Institute of Biology III, Faculty of Biology, University of Freiburg, 79104 Freiburg, Germany



dx.doi.org/10.17504/protocols.io.wj5fcq6

CyanoWorld



ABSTRACT

100x BG11 (500 ml)

74.79 g NaNO₃ $3.75 \text{ g MgSO}_4 \times 7H_2O$ 1.8 g CaCl₂ × 2H₂O 0.30 g Citric Acid 0.28 ml 0,5 M Na₂EDTA (pH 8)

Add Aqua bidest (ddH20) to a final volume of 500 ml. Aliquot (5x100 ml) and autoclave. Store at 4°C.

Trace-Metal-Mix (500 ml)

1.43 g H₃BO₃ 900 mg $MnCl_2 \times 4H_2O$ $110 \text{ mg ZnSO}_4 \times 7H_2O$ $195 \text{ mg Na}_2\text{MoO}_4 \times 2\text{H}_2\text{O}$ $39.5 \text{ mg CuSO}_4 \times 5H_2O$ $24.7 \text{ mg Co}(NO_3)_2 \times 6H_2O$

Put some ddH₂O into the bottle; dissovle the salts seperately and fill up to 500 ml. Sterile filter and aliquot the solution (10 x 50 ml). Store at 4°C.

2x BG11 (1 L)

20 ml 100x BG11 $2 \text{ ml } \text{K}_2 \text{HPO}_4 \times 3 \text{H}_2 \text{O (30 mg/ml)}$ 2 ml Na₂CO₃ (20 mg/ml) 2 ml Trace-Metal-Mix

Add 900 ml deionized water, afterwards add the solutions, fill up to 1 L with deionized water and autoclave. Store at RT. Before use add 2 ml of Fe-Ammonium-Citrate (6 mg/ml).

1xBG11 (1 L)

10 ml 100xBG11 10 ml TES pH 8 1 ml K₂HPO₄ (30 mg/ml) 1 ml Na₂CO₃ (20 mg/ml) 1 ml Trace-Metal-Mix

Add 900 ml deionized water, afterwards add the solutions, fill up to 1 L with deionized water and autoclave. Store at RT.

Before use add 1 ml of Fe-Ammonium-Citrate (6 mg/ml).

This is an open access protocol distributed under the terms of the Creative Commons Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited