



DNA extraction from whatman filter papers

Xiang Guo¹, Xiaoqing Zhang¹, Xiaohong Zhou¹

¹Department of Pathogen Biology, Key Laboratory of Prevention and Control for Emerging Infectious Diseases of Guangdong Higher Institutes, Guangdong Provincial Key Laboratory of Tropical Disease Research, School of Public Health, Southern Medical University

[dx.doi.org/10.17504/protocols.io.vjme4k6](https://doi.org/10.17504/protocols.io.vjme4k6)

Prof. Zhou Research Group



Xiang Guo

Department of Pathogen Biology...



ABSTRACT

This protocols is for RNA extraction from Whatman filter paper.

PROTOCOL STATUS


Working



We use this protocol in our group and it is working

STEPS MATERIALS



NAME	CATALOG #	VENDOR
Chelex 100	C7901-100G	Sigma Aldrich
DNase/RNase free distilled water	10977023	Thermo Fisher Scientific

Lyse RBC

- 1 Filter paper disc containing blood cut and placed into 1.5 mL EP tubes.
- 2 Add  1 ml RNase-free water and vortex to mix.

DNase/RNase free distilled water
 by [Thermo Fisher Scientific](#)
 Catalog #: 10977023
- 3 Incubate at  25 °C Room temperature for  00:15:00 to lyse RBC.

Precipitate the DNA and separate phases

- 4 Centrifuge at 20,000 × g for  00:05:00, discard the supernatant with a micropipettor.
 [go to step #2](#) until the DNA precipitate forms a white or pink gel-like pellet at the bottom of the tube.

5 Resuspend the pellet in 10% Chelex-100 solution, vortex the sample briefly for .



Chelex 100
by [Sigma Aldrich](#)
Catalog #: C7901-100G

6 Incubate the sample at for .

7 Centrifuge for at 20,000 × g.

8 Carefully transfer the supernatant to a clean 1.5 mL EP tube with a micropipettor.
Proceed to downstream applications, or store the DNA at .

Determine the DNA yield

9 Measure absorbance at 230nm, 260nm, and 280nm. Calculate the A260/A280 and A230/A260 ratio.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited