

# RNAqueous with TURBO DNA-free and SUPERase-In

Daniel Richter

## Abstract

Daniel Richter, Nov 29, 2011

based on RNAqueous May 29, 2008 protocol revision C, TURBO DNA-free June 9, 2009 protocol 1907M revision F

**Citation:** Daniel Richter RNAqueous with TURBO DNA-free and SUPERase-In. **protocols.io**

dx.doi.org/10.17504/protocols.io.iqzcdx6

**Published:** 04 Dec 2017

## Guidelines

To avoid possible RNA degradation, try to work quickly in all steps.

## Protocol

### Elution Solution Preparation

#### Step 1.

Place elution solution in a dry bath at 70-80°C

 **AMOUNT**

60 µl Additional info: per pellet

### Sample Lysate Preparation

#### Step 2.

Place pellet in a 50 mL conical tube on ice

### Sample Lysate Preparation

#### Step 3.

Add lysis buffer

 **AMOUNT**

700 µl Additional info:

### Sample Lysate Preparation

#### Step 4.

Pipette up and down 5 times to resuspend pellet

## Sample Lysate Preparation

### Step 5.

While still on ice, add additional lysis buffer

 [AMOUNT](#)

700 µl Additional info:

## Sample Lysate Preparation

### Step 6.

Pipette up and down 20 times to fully resuspend/lyse pellet

## Sample Lysate Preparation

### Step 7.

Transfer into two 1.5 mL Eppendorf tubes on ice

 [AMOUNT](#)

700 µl Additional info: each

## Sample Lysate Preparation

### Step 8.

Aliquot onto glass slide to check lysis from one Eppendorf tube

 [AMOUNT](#)

10 µl Additional info:

 [NOTES](#)

**Alyssa Alsante** 01 Jul 2017

Keep the slide inside a petri dish to ensure the noxious smell is contained

## Sample Lysate Preparation

### Step 9.

Immediately centrifuge both tubes at max speed in a microcentrifuge at 1°C, 5 min to pellet unlysed bacteria

## Sample Lysate Preparation

### Step 10.

During centrifugation, add cover slip to glass slide and quickly confirm lysis using a microscope

 [NOTES](#)

**Alyssa Alsante** 01 Jul 2017

If cells were lysed, transfer supernatant into two new 1.5 mL Eppendorf tubes; otherwise, continue with the same tubes

## Sample Lysate Preparation

## Step 11.

On ice, pass the lysate in each tube 5 times through a 25 gauge needle fitted to an RNase-free 1 mL syringe

### 📌 NOTES

**Alyssa Alsante** 01 Jul 2017

If the cells were not lysed in step 2, pass through the needle an additional time, and re-check lysis before continuing.

Filter Binding, Washing and Elution of RNA

## Step 12.

Add 64% Ethanol (at room temperature (RT)) per tube

### 📌 AMOUNT

700 µl Additional info:

Filter Binding, Washing and Elution of RNA

## Step 13.

Mix thoroughly by inverting the tubes several times

Filter Binding, Washing and Elution of RNA

## Step 14.

Transfer lysate/ethanol mixture from one of the Eppendorf tubes to a single filter cartridge in a collection tube

### 📌 AMOUNT

700 µl Additional info:

Filter Binding, Washing and Elution of RNA

## Step 15.

Centrifuge at 15,000 xg, RT, 30 sec or until the mixture is drawn through the filter

Filter Binding, Washing and Elution of RNA

## Step 16.

Discard the flow through

Filter Binding, Washing and Elution of RNA

## Step 17.

Repeat steps 14-16 three additional times until all of the lysate/ethanol mixture in both tubes has been drawn through the filter cartridge

Filter Binding, Washing and Elution of RNA

## Step 18.

Add wash solution #1

#### AMOUNT

700 µl Additional info:

Filter Binding, Washing and Elution of RNA

#### **Step 19.**

Centrifuge at 15,000 xg, RT, 30 sec

Filter Binding, Washing and Elution of RNA

#### **Step 20.**

Discard the flow through

Filter Binding, Washing and Elution of RNA

#### **Step 21.**

Add wash solution #2/3

#### AMOUNT

500 µl Additional info:

Filter Binding, Washing and Elution of RNA

#### **Step 22.**

Centrifuge at 15,000 x g, RT, 30 sec

Filter Binding, Washing and Elution of RNA

#### **Step 23.**

Discard the flow through

Filter Binding, Washing and Elution of RNA

#### **Step 24.**

Repeat wash steps 21-23

Filter Binding, Washing and Elution of RNA

#### **Step 25.**

Centrifuge an additional 30 sec at 15,000 x g

#### NOTES

**Alyssa Alsante** 03 Jul 2017

This step removes residual wash solution #2/3

Filter Binding, Washing and Elution of RNA

#### **Step 26.**

Place filter cartridge in a new collection tube

#### Filter Binding, Washing and Elution of RNA

##### Step 27.

Pipette pre-heated elution solution onto the filter

 **AMOUNT**

40 µl Additional info:

#### Filter Binding, Washing and Elution of RNA

##### Step 28.

Centrifuge at 15,000 xg, 30 sec

#### Filter Binding, Washing and Elution of RNA

##### Step 29.

Pipette an additional pre-heated elution solution onto the filter

 **AMOUNT**

10 µl Additional info:

#### Filter Binding, Washing and Elution of RNA

##### Step 30.

Centrifuge at 15,000 x g, 30 sec

 **NOTES**

**Alyssa Alsante** 03 Jul 2017

Final volume of eluate should be ~45 uL

#### TURBO DNase Treatment

##### Step 31.

Spec on NanoDrop and make a dilution at 1:5 (2 µl :8 µl) in a separate tube for PCR/Bioanalyzer

 **AMOUNT**

2 µl Additional info: aliquot

#### TURBO DNase Treatment

##### Step 32.

Take an aliquot of 1:5 dilution from previous step and dilute to 1 ng/µL in a separate tube (for Bioanalyzer)

 **NOTES**

**Alyssa Alsante** 03 Jul 2017

If concentration is >240 ng/uL (10 ug / 41.5 uL), dilute an aliquot to 240 ng/uL in 41.5 uL and save remainder at -80 C

## TURBO DNase Treatment

### Step 33.

Add SUPERase-In (at 20 U/μL, final concentration at 1 U/μL)

📄 AMOUNT

2.5 μl Additional info:

## TURBO DNase Treatment

### Step 34.

Add 10X TURBO DNase Buffer

📄 AMOUNT

1 μl Additional info:

## TURBO DNase Treatment

### Step 35.

Incubate at 37°C for 30 min

## TURBO DNase Treatment

### Step 36.

Vortex DNase Inactivation Reagent to resuspend

## TURBO DNase Treatment

### Step 37.

Add DNase Inactivation Reagent

📄 AMOUNT

5 μl Additional info:

## TURBO DNase Treatment

### Step 38.

Incubate 5 min at 22-26°C, flicking to mix 2-3 times during incubation.

## TURBO DNase Treatment

### Step 39.

Centrifuge at max speed for 1.5 min at RT

## TURBO DNase Treatment

### Step 40.

Transfer RNA to a fresh tube

📄 AMOUNT

50 μl Additional info:

📌 NOTES

To avoid transferring beads, pipette out at an angle and rotate tube such that beads are on the bottom.

#### TURBO DNase Treatment

##### **Step 41.**

Repeat centrifugation/transfer (steps 39-40) 2-3 times until no beads are transferred

#### TURBO DNase Treatment

##### **Step 42.**

Spec on NanoDrop and dilute 1  $\mu$ L aliquot of RNA to 1 ng/ $\mu$ L in a separate tube (for Bioanalyzer)