

# Measuring Droplet Volume in Home-Made Microfluidic Devices

Steve McCarroll's lab

## Abstract

Drop-seq is a technology we developed to enable biologists to analyze RNA expression genome-wide in thousands of individual cells at once. We first described this in [a 2015 paper in Cell](#). Though commercial implementations of droplet-based single-cell RNA-seq also now exist, we have made Drop-seq open-source and want to make sure that any lab can build their own system. The materials for constructing a Drop-seq setup in one's own lab can be obtained for about \$6,000. The reagents for performing Drop-seq cost about 6 cents per cell.

This is a supplemental protocol of our [Drop-Seq Protocol](#) for measuring droplet volume in home-made microfluidic devices.

**Citation:** Steve McCarroll's lab Measuring Droplet Volume in Home-Made Microfluidic Devices. **protocols.io**  
dx.doi.org/10.17504/protocols.io.mj5c4q6

**Published:** 23 Mar 2018

## Guidelines

### Videos and FAQs

These [tutorials, images, and diagrams](#) may be helpful in building your own Drop-seq setup and doing Drop-seq experiments in your lab.

We have also created a YouTube channel with a variety of [videos](#) to help scientists through the steps that most benefit from watching.

This [FAQ](#) provides also provides supplementary information.

## Before start

To measure droplet volume, purchase some durable, monodisperse polystyrene beads with a hydrophilic coating (e.x. 10-micron carboxylated polystyrene beads from Bangs Labs, product #PC06N-11355. It can be helpful to use fluorescent beads to be sure you can identify them in droplets. Bangs cells these under product # FC06F-10163).

## Protocol

### Step 1.

Wash and resuspend beads in Drop-seq lysis buffer at a concentration of 1000 beads per microliter.

### Step 2.

Draw the beads into a syringe with a magnetic mixer (as you would with the standard barcoded beads) and load into a syringe pump.

### Step 3.

Load the syringe pump intended for cells with regular PBS.

#### 📌 NOTES

**Anita Bröllochs** 11 Jan 2018

Since we are co-flowing beads with PBS, we estimate that the concentration of beads in the droplet fluid will be 500 beads per microliter

### Step 4.

Connect all tubing to the appropriate channels in the microfluidic device, and generate droplets.

### Step 5.

For a given number of droplets, count the number of beads inside. You should count the beads inside several hundred droplets to make sure that you have a statistically sound estimate.

### Step 6.

Divide the total number of beads counted inside droplets by the number of droplets you counted. This is your **droplet occupancy**.

### Step 7.

Here is how to calculate droplet volume:

*Droplet volume = (droplet occupancy) / (500 beads per microliter ) = # microliters per droplet.*