

# Nile Red Staining of *Drosophila* Larval Tissues

Elizabeth Allen

## Abstract

This protocol is used to stain late larval *Drosophila* lipid droplets in fat bodies and intestines with Nile Red, which emits fluorescence in the 552/636 nm range.

**Citation:** Elizabeth Allen Nile Red Staining of *Drosophila* Larval Tissues. **protocols.io**

dx.doi.org/10.17504/protocols.io.fdnbi5e

**Published:** 22 Jul 2016

## Guidelines

In brief, dissect animals in ice cold phosphate-buffered saline (PBS). Keep tissues in PBS on ice while obtaining your desired sample size. Fix tissues in 4% PFA (diluted in PBS with 0.1% Triton X-100), for 20 minutes, wash 3x in PBS, and stain tissues light-protected at room temperature for 1 hour.

\* Do not use any serums for this protocol because Nile Red will instead be drawn away from your tissues and into serum.

## Before start

1. Prepare Nile Red (Sigma-Aldrich) in acetone (1000 ug/mL).
2. Dissect tissues in ice cold PBS.
3. Fix tissues in 4% paraformaldehyde (formalin) diluted in PBS with 0.1% Triton X-100.
4. Wash tissues 3x in PBS.

## Materials

- ✓ 1X PBS (Phosphate-buffered saline ) by Contributed by users

- ✓ Triton-X100 by Contributed by users
- ✓ Paraformaldehyde Powder (PFA) [P6148](#) by Contributed by users
- Nile Red [N3013 SIGMA](#) by [Sigma Aldrich](#)
- Acetone solution [48358 SUPELCO](#) by [Sigma Aldrich](#)

## Protocol

### Step 1.

Dissect tissues in ice cold PBS, keeping samples on ice until required sample size is obtained.



#### REAGENTS

- ✓ 1X PBS (Phosphate-buffered saline ) by Contributed by users

### Step 2.

Fix tissues in 4% PFA for 20 minutes.



#### DURATION

00:20:00

### Step 3.

Wash tissues 3x in PBS.



#### REAGENTS

- ✓ 1X PBS (Phosphate-buffered saline ) by Contributed by users

### Step 4.

Stain samples in Nile Red at 0.5 ug/mL diluted in PBS for 1 hour.



#### DURATION

01:00:00



#### PROTOCOL

#### . [Nile Red prepared in acetone](#)

CONTACT: [Elizabeth Allen](#)

#### Step 4.1.

Prepare a concentrated working solution of Nile Red in acetone at 1000 ug/mL.



#### REAGENTS

Nile Red [N3013 SIGMA](#) by [Sigma Aldrich](#)

#### Step 4.2.

Store concentrated Nile Red solution at 4°C in the dark for up to 3 months.

### Step 4.3.

Use Nile Red/acetone concentrate diluted in PBS at a concentration of 0.5 ug/mL.

\*adjust duration of staining according to the tissue type, and stain at room temperature in the dark.

### Step 5.

Wash tissues 3x in PBS.



#### REAGENTS

✓ 1X PBS (Phosphate-buffered saline ) by Contributed by users

### Step 6.

Carefully replace the PBS with mounting medium before transferring samples to slides and imaging.

\* Image immediately, or temporarily store slides at 4°C.



#### PROTOCOL

#### . [Mounting Media for Immunohistochemistry - Drosophila](#)

CONTACT: [Sonia Hall](#)

#### Step 6.1.

90% glycerol

10% 1M Tris-base pH 8.0

0.5% n-propyl-gallate