

Western blotting in *Chlamydomonas reinhardtii*

João Vitor Molino

Abstract

This protocols describe the steps to perform a western blot in *Chlamydomonas reinhardtii* cell lysate and supernatant samples.

Citation: João Vitor Molino Western blotting in *Chlamydomonas reinhardtii*. **protocols.io**

[dx.doi.org/10.17504/protocols.io.kfpctmn](https://doi.org/10.17504/protocols.io.kfpctmn)

Published: 01 Nov 2017

Guidelines

Prepare cell material by cultivating in liquid media.

Before start

- Check antibody dilutions
- Check buffers disponibility
- Prepare 5% milk solution

Protocol

Sample preparation - Supernatant

Step 1.

1. Centrifuge algae culture at 2000xg for 10 min
2. Recover supernatant

Sample preparation - Lysate

Step 2.

1. Centrifuge algae culture at 2000xg for 10 min
2. Remove supernatant, and re-suspend cells in lysis buffer (50 mM Tris·HCL (pH 8.0), 0.1% Triton X-100), concentrating cells 100-fold.

Sonication

Step 3.

1. Sonicate using appropriate sonication tip. Duty cycle: 0.5 | Amplitude: 20% | Duration: 30 s

2. Centrifuge for 15 min at 20000xg to remove cells debris and recover soluble proteins
3. Quantificate soluble protein

Gel electrophoresis | SDS-PAGE

Step 4.

1. Load 30 µg of total soluble protein (TSP) per lane in a 12% SDS-PAGE
2. Transfer proteins to a nitrocellulose membrane
3. Block the membrane with 5% milk solution
4. Probe the desired protein with the specific antibody, diluted in 5% milk solution
5. Wash 2 times with TBST (0.2 M Tris, 1.37 M NaCl, 0.1% Tween-20, pH 7.6)
6. Add secondary antibody if required

AMOUNT

30 µg Additional info: total soluble protein per lane