



May 17,
2019

Working

Yale - Triglycerides [↗](#)

John Stack¹, Gary Cline¹

¹Yale University

[dx.doi.org/10.17504/protocols.io.y4bfysn](https://doi.org/10.17504/protocols.io.y4bfysn)

Mouse Metabolic Phenotyping Centers
Tech. support email: info@mmpc.org

Lili Liang

ABSTRACT

Summary:

Procedure used to determine the concentration of triglycerides in blood, serum, and plasma. Triglycerides are determined by coupling lipase, glucokinase, glycerol phosphate oxidase, and peroxidase to form a quinonemine dye which is measured at 500 nm.

EXTERNAL LINK

<https://mmpc.org/shared/document.aspx?id=213&docType=Protocol>

MATERIALS

NAME	CATALOG #	VENDOR
Triglyceride Reagent and Activator	R85457	Cliniqa
Multi Analyte Calibrator	R60010	Prolabs(cliniqa)
Assayed Control Serum 1	R83082	Prolabs(cliniqa)
Assayed Control Serum 2	R83083	Prolabs(cliniqa)

MATERIALS TEXT

Reagent Preparation:

Triglyceride Reagent: Add 40mL of Triglyceride Activator to the Reagent bottle.

Multi Analyte Calibrator: Add the appropriate amount of water (6.5mL) to the chemical control bottle. Invert to mix, allowing 15 minutes for the reagent to settle.

Assayed Control Serum 1: Add the appropriate amount of water (6.5mL) to the chemical control bottle. Invert to mix, allowing 15 minutes for the reagent to settle.

Assayed Control Serum 2: Add the appropriate amount of water (6.5mL) to the chemical control bottle. Invert to mix, allowing 15 minutes for the reagent to settle.

BEFORE STARTING

Analysis by automated system Cobas Mira Plus

1 Calibrate Cobas for Triglyceride analysis by running a multi-analyte calibrator and two control serum.

2 Sample handling as performed by Cobas Mira Plus.

- a) Pipette 4 μL of sample into cuvette.
- b) Add 275 μL of Triglyceride liquid reagent.
- c) Incubate at 37°C for 10 minutes.
- d) Absorbance is measured at 500nm.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited