# Cell lysis, detergent-free

#### **Queenie Chan**

#### **Abstract**

Detergents are generally not compatible with mass spectrometers, so this is a detergent-free method of cell lysis that is compatible with mass spectrometry. Since this protocol does not have a precipitation step, it saves time and minimizes sample loss as well.

Citation: Queenie Chan Cell lysis, detergent-free. protocols.io

dx.doi.org/10.17504/protocols.io.k79czr6

Published: 02 Mar 2018

#### **Guidelines**

Mix well after each step by vortex for 2 seconds.

## **Before start**

Prepare ice bucket.

Prepare 90°C heat block.

#### **Protocol**

#### Cell sample

#### Step 1.

Fresh or previously frozen, store on ice. Original protocol uses  $5x10^7$  cells per sample.

# Resuspend cells

#### Step 2.

**■** AMOUNT

1.5 ml Additional info: ice cold Milli-Q water

Add trifluoroethanol (TFE)

Step 3.

**■** AMOUNT

1.5 ml Additional info: Trifluoroethanol

**A** SAFETY INFORMATION

This step should be done in the fume hood.

NOTES

Queenie Chan 11 Dec 2017

1:1 water-TFE acts as a hypotonic aqueous buffer to lyse cells, eliminating the need for detergent. TFE helps protein solubility and denaturation; it readily evaporates, so removing it is easy.

Queenie Chan 11 Dec 2017

TFE is located in NCE 435 - flammable cabinet

Queenie Chan 11 Dec 2017

TFE evaporates fast, so work quickly.

Sit on ice

Step 4.

Vortex

Step 5.

Sonicate

Step 6.

#### NOTES

Queenie Chan 11 Dec 2017

Original protcol calls for 2 mins at 20% duty cycle, output control 3 using a ultrasonic probe. The time required to do this for large number of samples is not realistic, and therefore not used here.

Queenie Chan 11 Dec 2017

While this step is happening, turn on the heat block and set it at 95oC

#### Adjust pH

Step 7.

**AMOUNT** 

200 µl Additional info: Tris pH 8.5, 100mM

#### Remove debris

Step 8.

Spindown debris at top speed on a benchtop centrifuge for 10 minutes, then keep the supernatant.

#### Measure protein concentration

#### Step 9.

Use the Coomassie reagent to measure protein concentration. Aliquot the desired amount to proceed with the remainder of this experiment.

#### NOTES

Queenie Chan 11 Dec 2017

Freezing unused protein is not recommended since there is a high chance of precipitation upon thawing.

#### Reducing agent

#### **Step 10.**

Add TCEP (tris(2-carboxyethyl)phosphine to a final concentration of 10mM

#### **P** NOTES

Queenie Chan 11 Dec 2017

TCEP is located at NCE 436 -20C freezer door

#### Alkylating agent

#### **Step 11.**

Add 2-chloroacetamide to a final concentration of 40mM

#### NOTES

Queenie Chan 11 Dec 2017

2-chloroacetamide is located in NCE 438 - chemical shelves

#### Preparing the lysate for digestion

#### Step 12.

**■** AMOUNT

15 ml Additional info: 50mM ammonium bicarbonate

#### Digestion with LysC

**Step 13.** 

Add 1ug LysC to every 50ug protein

#### **▮** TEMPERATURE

37 °C Additional info: incubator (not heat block)

#### Digestion with trypsin

#### Step 14.

Add 1ug trypsin to every 50 ug protein. Allow reaction to occur overnight.

# **▮** TEMPERATURE

37 °C Additional info: incubator (not heat block)

## STAGE tip

Step 15.

Proceed with desalting as normal