

Western blot

Laura Ruiz Remolina

Abstract

Citation: Laura Ruiz Remolina Western blot. **protocols.io**

dx.doi.org/10.17504/protocols.io.j8bcrsn

Published: 07 Oct 2017

Guidelines

Avoid repeated freeze-thaw cycles

Before start

Prepare all the reagents and mind the pH of each one

Materials

- Phosphatase inhibitor PhosSTOP 04 906 837 by [Roche](#)
- Proteases inhibitor Complete 11 697 498 by [Roche](#)
- DC protein assay 5000112 by [BIO-RAD](#)
- BSA A7906 by [Sigma Aldrich](#)
- Acrlamide 161-0107 by [BIO-RAD](#)
- 2-mercaptoethanol M-6250 by [Sigma-aldrich](#)
- APS A-3678 by [Sigma Aldrich](#)
- Bisacrilamide 161-0201 by [BIO-RAD](#)
- H2O2 8597 by [Merck Millipore](#)
- Iodophenol I1020-1 by [Sigma-aldrich](#)
- Glicerol 4095 by [Merck Millipore](#)
- Glycine 161-0724 by [BIO-RAD](#)
- ✓ HCl 20248.295 by Contributed by users
- Luminol A-8511 by [Sigma-aldrich](#)
- ✓ PVDF membranes #66543 by Contributed by users
- ✓ Methanol 141091.1214 by Contributed by users

NaCl S-3014 by [Sigma-aldrich](#)

TEMED (Tetramethyl-ethulenediamine) T9281 by [Sigma-aldrich](#)

SDS 161-0302 by [BIO-RAD](#)

Tween-20 P-7949 by [Sigma-aldrich](#)

Protocol

- Preparation of lysate from tissue—Remove a small volume of lysate to perform a protein quantification assay (with DC protein assay, BIO-RAD). Determine the protein concentration for each tissue lysate.

Step 1.

Load equal amounts of protein into the wells of the SDS -PAGE gel, along with molecular weight marker—60 ug total proteins per pore—Run the gel for 90min at 100 V—Gel percentage —separation gel 10%—spacer gel 5%.

Step 2.

Activate PVDF with methanol for 1 min and rinse with transfer buffer before preparing the stack—Run at 100 V for 90min.

Step 3.

- Wash the membrane in three washes of TBST—5 percent skimmed milk powder—5 min each

Step 4.

Block the membrane 4 hours in BSA 5%. Wash the membrane in 3 washes of TBST, 5 min each.

Step 5.

Wash 3 times, 5 minutes each

Step 6.

Incubate the membrane with the primary antibody diluted as the manufacturer's instructions.

Step 7.

Wash 3 times, 5 minutes each

Step 8.

Incubate the membrane with the recommended dilution of conjugated secondary antibody in blocking buffer at room temperature for 1 h

Step 9.

. Wash the membrane in 3 washes of TBST, 5 min each

Step 10.

Acquire image using normal image scanning methods for colorimetric detection.

Step 11.