



Dec 20, 2019

Dye-terminator DNA sequencing V.2

Diep R. Ganguly¹

¹The Australian National University

1 Works for me dx.doi.org/10.17504/protocols.io.baq5idy6

Pogson Genomics Group



Diep R. Ganguly
The Australian National University



ABSTRACT

This protocol (based on the BigDye® Terminator v3.1 Cycle Sequencing Kit) is for performing terminator cycling sequencing reactions for Sanger sequencing of amplified PCR products or plasmid DNA on the 3130X genetic analyser (Applied Biosystems).

□
BigDye Terminator
v3.1.pdf

MATERIALS

NAME	CATALOG #	VENDOR
Antarctic Phosphatase - 1,000 units	M0289S	New England Biolabs
96 well PCR Plate Non-skirted	MPS-499	Phenix Research
Nuclease-free water (e.g. MilliQ or HPLC grade water)		
primers		
EDTA		
10 mM dNTPs	10297-018	Life Technologies
Ethanol	100983	Merck Millipore
BigDye™ Terminator v3.1 Cycle Sequencing Kit	4337454	Thermo Fisher
Exonuclease I (E. coli)	M0293S	NEB
Hi-Di™ Formamide	4311320	Thermo Fisher Scientific

BEFORE STARTING

Optimize PCR cycling (if sequencing amplified PCR products) to ensure your reaction produces a single product. Perform gel excision or PCR clean-up with the potential inclusion of incubating with Antarctic phosphatase and Exonuclease 1 to dephosphorylate and degrade unincorporated dNTPs in PCR reactions to prepare templates for DNA sequencing.

Terminator cycling reaction

- 1 Perform sequencing reaction in PCR tubes (or 96-well plate) with BigDye Terminator cycling kit and forward or reverse primers.^{30m}

Component	Volume (µl)
2.5X Reaction Ready Mix	1
5X BigDye Sequencing buffer	2
20 µM F/R Primer	0.5
Template (plasmid or cleaned PCR product)	150 ng dsDNA
Nuclease-free water	to 10 µl

BigDye Terminator Cycling reaction

5x reaction buffer: 400 mM TRIS, 10 mM MgCl₂

- 2 Run the following thermal cycling protocol:
 1. 1 min at 96 °C
 2. 30-40 cycles: 96 °C for 10 seconds, 50 °C for 5 seconds, and 60 °C for 4 min.
 3. Hold at 4-12 °C.

4h

Purification 1h 30m

- 3 Transfer PCR reaction to eppendorf tube. To the reaction, add 2.5 µL of 125 mM EDTA (make sure it touches bottom of tube).
- 4 Add 30 µl of 100% ethanol, mix well (inversion).
- 5 Incubate at room temperature for 15 minutes.
- 6 Centrifuge at 4 °C at max speed for 30 minutes.
- 7 Discard supernatant and add 50 µl of ice-cold 70% ethanol.
- 8 Centrifuge at 4 °C at max speed for 5 minutes.
- 9 Discard supernatant and allow to air-dry in the dark for >15 minutes.

Prepare for sequencing

- 10 Resuspend the pellet (likely transparent) in 7.5 µL HiDi Formamide (add to any empty wells). Incubate at RT for 5 minutes then transfer to plate. Spin down briefly.
- 11 Incubate plate at 95 °C for 3 minutes (denature) then place immediately on ice. Spin down briefly.
- 12 Submit for sequencing on 3130X genetic analyser (Applied Biosystems). Keep samples on ice.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited