

Splenocyte Preperation

Rabia Khan

Abstract

Citation: Rabia Khan Splenocyte Preperation. **protocols.io**

dx.doi.org/10.17504/protocols.io.es9beh6

Published: 20 Oct 2016

Protocol

Step 1.

Harvest mouse spleen and prepare a single cell suspension. Use slides or a syring to push the spleen through a cell strainer.

Step 2.

Pellet the cells by centrifugation (350 x g); aspirate the supernatant.

Step 3.

Dilute the 10X Red Blood Cell Lysis Buffer to 1X working concentration with deionized water and resuspend the pellet in 5 ml of 1X Lysis Buffer.

Step 4.

Incubate on ice for 4-5 minutes with occasional shaking.

Step 5.

Stop the reaction by diluting the Lysis Buffer with 20-30 ml of 1X PBS.

Step 6.

Spin the cells (350 x g) and discard the supernatant.

Step 7.

Resuspend the pellet in the appropriate buffer

Step 8.

Count cells, adjust density, and proceed with cell staining procedures.