



May 14,
2019

Working

UC Davis - Total Cholesterol (TC) Protocol [↗](#)

Peter Havel¹

¹University of California, Davis

[dx.doi.org/10.17504/protocols.io.ygeftte](https://doi.org/10.17504/protocols.io.ygeftte)

Mouse Metabolic Phenotyping Centers
Tech. support email: info@mmpc.org

Lili Liang

ABSTRACT

Summary:

Cholesterol esters are enzymatically hydrolysed by cholesterol esterase to cholesterol and free fatty acids. Free cholesterol, including that originally present, is then oxidized by cholesterol oxidase to cholest-4-en-3-one and hydrogen peroxide. The hydrogen peroxide combines with HBA and 4-aminoantipyrine to form a chromophore (quinoneimine dye) which may be quantitated at 500-550nm.

EXTERNAL LINK

<https://mmpc.org/shared/document.aspx?id=92&docType=Protocol>

MATERIALS

NAME	CATALOG #	VENDOR	CAS NUMBER	RRID
Calibrator	TR43002	Fisher Diagnostics		
Reagents	TR13421	Fisher Diagnostics		
PBS				
Microplate				
Platereader				

MATERIALS TEXT

Note:

Thermo Fisher Scientific, [RRID:SCR_008452](#)

- 1 Add 5 µl of calibrator and sample to each well.

IMPORTANT: Make sure not to add any bubbles to the wells when dispensing reagents, this will interfere with reading in the platereader.

- 2 Add 300 µl of reagent to each well. Incubate at 37°C for 5 minutes. Read at 540 nm.

IMPORTANT: If samples are hemolyzed, pipet a blank well with 5 µl sample and 300 µl PBS

- 3 Subtract blank readings from final readings. The assay will be linear so the unknown samples can be calculated as $(\text{Sample Absorbance} \div \text{Calibrator Absorbance}) \times \text{Calibrator Concentration}$.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited