



Nov 17,
2019

Salty sample derivatization protocol for GC-MS

Forked from [Salty sample derivatization protocol for GC-MS \(BSTFA\)](#)

[Emilia Sogin](#)¹, [Erik Pukus](#)¹, [Nicole Dubilier](#)¹, [Manuel Liebeke](#)¹

¹Max Planck Institute for Marine Microbiology

1 *Works for me* [dx.doi.org/10.17504/protocols.io.nyxdfxn](https://doi.org/10.17504/protocols.io.nyxdfxn)

Mass Spectrometry at MPI-Bremen

[Emilia Sogin](#)
Max Planck Institute for Marine Microbiology

ABSTRACT

Metabolomic footprinting, or exometabolomics, seeks to describe the complex chemical space surrounding cells. However, metabolite detection is hampered by negative matrix effects, especially when samples originate from high salt-content environments. We propose a sample preparation method for the untargeted exometabolomic analysis of seawater by gas chromatography-mass spectrometry. By investigating the independent effects of salt and water on metabolite ionization, we demonstrate residual water locked within salt crystals inhibits compound detection. Our novel preparation protocol includes a combination of solvent drying and ultra-sonication to facilitate the removal of water molecules from the salt palette and extraction of small molecules into the derivatization reagents. We can now detect a broad range of sugars, sterols and amino and organic acids from as little as 0.5 mL of seawater in the μM range. To demonstrate the utility of our approach, we monitored the exo-metabolome of *Marinobacter adhaerens*, a common marine bacterium during its exponential growth phase. With our untargeted method, we show *M. adhaerens* preferentially metabolizes glutamic acid and proline before valine, isoleucine and leucine. Our novel untargeted footprinting approach advances the capacity to detect small compounds from high salt content samples by filling a gap in analytical chemistry for the GC-MS analysis of marine samples.

MATERIALS

NAME	CATALOG #	VENDOR
Methoxyamine		CS-Chromatographie-Service GmbH
N,O-Bis(trimethylsilyl)trifluoroacetamide		CS-Chromatographie-Service GmbH

Sample Drying


- 1 Dry 0.5 to 1 ml of a salty sample for 6 hr using an eppendorf Concentrator plus SpeedVac at 45 °C in V-AQ mode.

Critical step: the drying in the SpeedVac should take no longer than 6h, otherwise the samples could degrade. Samples can be kept in SpeedVac overnight after the run

- 2 After removal from the SpeedVac, resuspend the sample in 250 μL of anhydrous toluene
- 3 Ultrasonicate the sample for 10 min at 100% efficiency in a sonication bath

- 4 Remove solvent under a steady flow of N₂ gas.

Important note: Prior to inserting needles into tubes, make sure they are pre-cleaned with organic solvent. Regulate gas flow so only the surface of the solvent is disturbed. Needles should not touch the solvent.




3-D printed drying device with attached needles.

Pause Point

- 5 Samples can be stored at 4°C until derivatization.

Derivatization

- 6 Add 80 µL of Methoxyamine (MeOx; 20 mg/1 mL of pyridine) to dried sample.



The salt crystals should dissolve into the pyridine creating a milky white solution as pictured.

- 7 Ultrasonicate for 10 min using a room temperature bath

- 8 Incubate the sample at 37°C for 90 minutes and 1350 r.p.m. using a thermomixer
- 9 Remove pyridine solvent from sample leaving behind salt crystals under a steady flow of N₂ gas. Use a thermomixer to maintain constant temperature (25°C) and rotational speed (200 r.p.m.).
- 10 Once sample is dry, add 100 µL of N,O-Bis(trimethylsilyl)trifluoroacetamide (BSTFA) to salt pellet.
Note: Do not leave tubes open for extended periods of time to avoid added moisture from the air
- 11 Ultrasonicate for 10 min using an ultrasonication bath at room temperature
- 12 Vortex the sample to dissolve salt crystals into the pyridine solvent
- 13 Incubate the sample at 37 °C for 30 min and 1350 r.p.m. using a thermomixer
- 14 Ultrasonicate the derivatized sample for 10 min
- 15 Pellet salt by centrifuging for 2 min at 14,800 rpm
- 16 Transfer 100 µL of the clear supernatant into a GC-MS vial insert and place on autosampler for GC-MS analysis

Critical step: the samples should be run within 1 day after the preparation, otherwise they will undergo degradation and the MeOX will crystallize, making the injection impossible



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited